INSTRUCTIONS FOR USE

Implants, abutments and associated restorative and dental laboratory components. For detailed information on the specific product system you are using, please consult the Adin Dental Implant Systems Ltd. catalog.

DESCRIPTION

The Adin dental implant system includes a variety of implants manufactured from titanium. Dental implants are surgically inserted into the upper and/or lower jawbone and serve as a substitute or replacement tooth root providing a stable foundation for restorations.

INDICATIONS FOR USE

Adin's implants are intended for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations including cement retained, screw retained, or overdenture restorations, and terminal or intermediate abutment support for fixed bridgework.

DIRECTIONS FOR USE

The implantation procedure should be done under aseptic conditions with specifically designed sterile surgical instruments, an electrical surgical drilling system with internal and external irrigation is recommended for drilling the surgical site. Specific drilling sequences for placement of implants should be followed.

WARNINGS

Improper technique can contribute to implant failure and / or loss of bone. Adin's dental implants are intended for use only in the indicated applications. Dental implants must not be altered in any way. The use of electro-surgical instruments or lasers around metallic implants and their abutments is not recommended due to the risk of electric and / or heat conductivity. Implant mobility, bone loss, or chronic infection may indicate implant failure. If the implant becomes contaminated by the patient's body fluids in any way, the implant cannot be used in any other patient.

PRECAUTIONS

Used uncleaned and unsterile devices may pose a biological hazard due to tissue contamination. To prevent risks associated with such hazards, dispose used devices in accordance with applicable local laws and regulations or according to institutional protocol. The surgical techniques required to place endosseous dental implants require specialized and complex procedures.

Formal training for placement of implants is recommended. Important: determine local anatomy and suitability at the available bone for implant placement. Adequate radiographs direct palpation and visual inspection of the implant site are necessary prior to treatment, planning and use of Adin implants.

ADVERSE REACTIONS

Some of the complications that can occur include: infection, bone loss, patient discomfort, implant mobility, local soft tissue degeneration, and unfavorable implant placement or alignment. Additional information and steps to be taken can be found in the User Manual.

CONTRAINDICATIONS

General contraindications associated with elective surgery should be observed. Important: determine local anatomy and suitability at the available bone for implant placement. Adequate radiographs direct palpation and visual inspection of the implant site are necessary prior to treatment, planning and use of Adin implants. Possible contraindications: chronic bleeding problems, psychological impairment, metabolic bone or connective tissue diseases, treatment with corticosteroids, certain cardiac and vascular diseases, tobacco usage, diabetes (uncontrolled), treatment with chemotherapeutic agents, chronic renal disease, poor patient oral hygiene, bruising alcoholism. Temporary contraindications: systemic infection, local oral or respiratory infection. Anatomical or pathological contraindications: insufficient alveolar bone width and height to surround the implant with at least one millimeter of bone, both buccally and lingually to the most superior aspect of the implant body, inadequate bone height where proper implant placement would encroach on the mandibular canal, malignancies. Unfavorable implant placement or alignment. Additional information and steps to be taken can be found in the User Manual.

HOW SUPPLIED

This package contains one dental implant. A cover screw is included with all two piece implants.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician. Adin dental implant has not been evaluated for safety and compatibility in the MR environment. Adin dental implant has not been tested for heating or migration in the MR environment. The maximum insertion Torque is 50 Ncm.

Do not reuse implants, cover screws, temporary abutments and abutments. Re-use of these items, will lead to an increased risk for product failure as functionality cannot be guaranteed if these products are reused. In addition, there is an increased risk of contamination.

EXPLANATION OF PICTOGRAMS

	Sterilized using irradiation		Manufacturer
	Batch number		Use by
	Catalogue number		CAUTION: See instructions for use
	Do not use if package is damaged		Do not re-use
	Do not resterilize		Date of manufacture
	Authorised representative in the European Community		ONLY

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■ Adin Dental Implant Systems LTD.
MedNet EC-REP GmbH