

CloseFit™

Regular & Wide Implants

CloseFit™


Regular & Wide Implants

Our strong CloseFit™ Conical Connection is designed for optimal seating while reducing and virtually eliminating micromovements.

Appropriate for immediate function/loading capability in all bone types.

CloseFit™ implants are designed to be used in all indications and in all regions of the lower and upper jaws.

The back-tapered coronal portion is designed to reduce pressure and enable healthy cortical bone around the implant collar .



CloseFit™ RP

3.5mmD
Regular Platform

* Built-In Platform
Switching

* Minimally
invasive surgery

* Osseofix™

L-8mm | L-10 mm | L-11.5 mm
L-13 mm | L-16 mm | L-18 mm


CloseFit™ WP

4.3/5.0mmD
Wide Platform


* Back
tapered collar

* Micro Threads:
Engineered to enhance bone
regeneration, prevent bone
loss & increase bone implant
contact

* Bone - condensing property
* Self - drilling




*The features related to WP & RP Implants




RP CloseFit™

Regular Platform



Recommended insertion torque up to 50 Ncm

CloseFit™ RP	Ø	L	D3	Catalog No.
	3.5	8	2.6	ISPX0835
		10	2.55	ISPX1035
		11.5	2.5	ISPX1135
		13	2.5	ISPX1335
		15	2.5	ISPX1535
		18	2.4	ISPX1835

Radiographic Template RP CloseFit™ | ASD0612


Drill Sequence

Recommended drilling speed up to 1200 rpm

Implant Diameter

ø3.5

ø3.5



You may use Tri-Step drill instead of ø2, ø2.8 and ø3.2 drills in sequence.


Cover Screw & Healing Abutments Recommended Torque 15(Ncm)

RP CloseFit™ Conical Cover Screw



RP0073

RP CloseFit™ Slim Healing Abutment



RP0046, GH2/L6.8


RP0047, GH3/L7.8

RP0048, GH4/L8.8

RP0049, GH5/L9.8

RP0050, GH6/L10.8

RP CloseFit™ Healing Abutment



RP0022, GH2/L6.8

RP0023, GH3/L7.8

RP0024, GH4/L8.8

RP0025, GH5/L9.8

RP0026, GH6/L10.8

RP CloseFit™ Conical Healing Abutment



RP0090, GH2/L7.2

RP0091, GH3/L8.2

RP0092, GH4/L9.2

RP0093, GH5/L10.2

RP0094, GH6/L11.2




RP™ Impression

Recommended Torque up to 15(Ncm)





Our components for taking impressions for all restoration types can be successfully utilized for digital or conventional protocols either implant-level or abutment level procedures.

Our broad range of solutions and flexible designs enable easy and accurate restorative planning, whether you're dealing with a single tooth, multiple teeth, bridge or full mouth rehabilitation.




Transfers		L	∅	D1	H	GH	RH	Catalog No.
	RP Slim Closed Tray Transfer - Short	11.2	3.6	2.9	9.0	-	-	RP5109
	RP Slim Closed Tray Transfer - Long	15.2	3.6	2.9	13.0	-	-	RP5113
	RP Closed Tray Transfer - Short	11.20	4.25	2.90	9.0	-	-	RP5209
	RP Closed Tray Transfer - Long	15.20	4.25	2.90	13.0	-	-	RP5213
	RP Slim Open Tray Transfer - Short	14.20	5.00	2.90	10.0	-	-	RP5610
	RP Slim Open Tray Transfer - Long	18.20	5.00	2.90	14.0	-	-	RP5614
	RP Open Tray Transfer - Short	14.20	5.00	2.90	10.0	-	-	RP5510
	RP Open Tray Transfer - Long	18.20	5.00	2.90	14.0	-	-	RP5514
	RP Scan Body	9.2	-	2.9	7.0	-	-	RP0088
		12.2	-	2.9	10.0	-	-	RP0089

Analog / Digital Analog		L	∅	D1	H	GH	RH	Catalog No.
	RP implant analog	12.60	3.55	2.90	-	-	-	RP0007
	RP Digital Model Analog 15(Ncm)	10.00	3.55	2.90	-	-	-	RP1000



Cement Retained Abutments

Recommended Torque 35(Ncm)





Abutments	L	Ø	D1	D2	H	GH	RH	Catalog No.
 RP Straight Slim Abutment No End Line	11.00	4.00	2.90	2.95	8.80	2.80	6.00	RP0020
 RP Esthetic Straight Abutment	9.70	4.40	2.90	3.00	7.50	1.00	6.50	RP0012
	10.70	4.40	2.90	3.00	8.50	2.00	6.50	RP0011
	11.70	4.40	2.90	3.00	9.50	3.00	6.50	RP0010
 RP Esthetic Angled Abutment 15°	9.70	4.40	2.90	2.34	7.50	1.00	6.50	RP0013
	10.70	4.40	2.90	2.34	8.50	2.00	6.50	RP0014
	11.70	4.40	2.90	2.34	9.70	3.00	6.50	RP0015
 RP Engaging Temporary Abutment	12.65	3.70	2.90	3.40	10.45	1.00	9.45	RP0043
 RP Non Engaging Temporary Abutment	12.65	3.70	2.90	3.40	10.45	1.00	9.45	RP0044
 RP Engaging Plastic Cylinder Abutment (Casting)	12.70	3.70	2.90	3.00	10.50	1.00	9.50	RP0045
 RP Engaging Titanium/Plastic Abutment (Casting)	13.70	4.70	2.90	4.70	11.50	1.50	10.00	RP3095
 RP Non-Engaging Titanium/Plastic Abutment (Casting)	13.70	4.70	2.90	4.70	11.50	1.36	10.00	RP3097



TMA™ - Multi Unit

Recommended Torque 35(Ncm)










Our TMA™ (Trans Mucosal Abutment) system is indicated for multiple-unit, screw-retained restorations and may be used in combination with an implant level framework design.

The TMA™ system is used to elevate the restorative seating platform when restoration at implant level is not indicated or practical due to the depth or angle of the implant.

The TMA™ System is designed to allow better prosthetics access when using implants with divergent angles accommodate full and partial edentulous arches, especially when using angled implants and all-on-x technique for full arch restoration.

The TMA™ abutment system is available in straight or angled (17°, 30°), engaging the Conical Connection for indexing, with a selection of collar heights. Abutments are delivered with handles for easy insertion and to assist with seating.

TMA		L	Ø	D1	H	GH	RH	Catalog No.
	RP Straight Trans Mucosal Abutment	4.20	4.90	2.90	2.00	1.00	2.10	RP0032
		5.20	4.90	2.90	3.00	2.00	2.10	RP0033
		6.20	4.90	2.90	4.00	3.00	2.10	RP0034
		7.20	4.90	2.90	5.00	4.00	2.10	RP3729
		8.20	4.90	2.90	6.00	5.00	2.10	RP3730
	RP Angled Trans Mucosal Abutment 17° 2.5mm	6.62	4.90	2.90	-	2.58	2.15	RP3734
	RP Angled Trans Mucosal Abutment 17° 3.5mm	7.65	4.90	2.90	-	3.68	2.15	RP3737
	RP Angled Trans Mucosal Abutment 30° 3mm	6.97	4.90	2.90	-	3.65	2.15	RP3735
	RP Angled Trans Mucosal Abutment 30° 4mm°	8.00	4.90	2.90	-	4.00	2.15	RP3738



TMA components		L	∅	D1	D2	H	GH	RH	Catalog No.
	TMA Healing Cap	5.5	4.9	-	3.5	-	-	-	RS5005
		7	4.9	-	3.5	-	-	-	RS5003
	TMA Impression Coping Open Tray	11.0	4.9	-	3.1	-	-	-	RS5006
	TMA Closed Tray Transfer	8.0	4.9	-	4.0	-	-	-	RS5011
	TMA Scan Body	6	4.9	4.9	3.8	-	-	-	TM0006
		9	4.9	4.9	3.8	-	-	-	TM0009
	TMA Titanium Temporary Cylinder	12.0	4.9	-	3.5	-	1.5	10.5	RS4900
	TMA Cementing Cone	4.9	4.9	-	2.9	-	1.5	3.4	TM0010
	TMA Plastic Casting Sleeve w/ Prosthetic Screw	12.0	4.9	-	3.3	-	1.5	10.5	RS5001
	TMA Abutment Replica	16.8	4.9	-	1.7	-	-	-	RS5004
	TMA Digital Model Replica Recommended torque 15(Ncm)	12.0	4.9	-	3.8	-	-	-	RS1001







Overdenture Restoration

Recommended Torque 35(Ncm)



Ball Attachment		L	Ø	D1	D2	H	GH	RH	Catalog No.
	RP Ball Attachment	8.86	3.10	2.90	2.50	4.00	0.50	3.52	RP0051
		9.32	3.10	2.90	2.50	4.50	1.00	3.52	RP0052
		10.32	3.10	2.90	2.50	5.50	2.00	3.52	RP0053
		11.32	3.10	2.90	2.50	6.50	3.00	3.52	RP0054
		12.32	3.10	2.90	2.50	7.50	4.00	3.52	RP0055
		13.32	3.10	2.90	2.50	8.50	5.00	3.52	RP0056


Ball Attachment caps		L	Ø	D1	D2	H	GH	RH	Catalog No.
	Plastic Ball Cap White (Nylon)	2.60	4.10	-	3.60	-	-	-	RS2660
	Plastic Ball Cap Pink Soft (Nylon)	2.60	4.10	-	3.60	-	-	-	RS2662
	Titanium Ball Cap	2.90	4.70	-	4.70	-	-	-	RS2670
	Stainless Steel Ball Cap	2.90	4.70	-	4.70	-	-	-	RS2675SS

Ball Attachment		L	Ø	D1	D2	H	GH	RH	Catalog No.
	Ball Attachment Abutment Replica	16.60	4.00	-	2.5	-	0.80	2.7	RS5740
















Screw - Retained

Recommended Torque 35(Ncm)



Flat Connection	L	Ø	D1	D2	H	GH	RH	Catalog No.
	6.80	4.50	2.90	-	2.04	2.04	-	RP0030
RP Flat Connection Abutment	8.80	4.50	2.90	-	4.04	4.04	-	RP0031
	11.5	4.5	-	3.9	11.1	0.4	10.7	FC5515
	12.6	4.5	-	3.4	-	-	-	FC5737
	13.6	4.6	-	3.0	13.6	0.1	13.5	FC6011
	5.0	2.5	-	-	-	-	-	FC6015
	3.8	4.5	-	3.1	3.4	0.6	2.8	FC6013
	17.0	2.5	-	-	-	-	-	FC6017
	3.8	4.5	-	3.2	3.4	0.6	2.8	FC6018
	11.0	4.5	-	4.2	-	-	10.2	FC5512

RP CloseFit™ Screws		L	Ø	Catalog No.
	RP Abutment Screw	7.20	2.20	RP0001
	TMA Prosthetic Screw 15(Ncm)	4.1	2.1	RS3404
	RP Retaining Screw for TMA 1mm	8.00	3.10	RP0039
	RP Retaining Screw for TMA 2mm	9.00	3.10	RP0040
	RP Retaining Screw for TMA 3mm	10.00	3.10	RP0041
	RP Retaining Screw for TMA 4mm	11.00	3.10	RP3409
	RP Retaining Screw for TMA 5mm	12.00	3.10	RP3412
	RP Retaining Screw for Angled TMA	7.00	2.35	RP3403
	RP Open Tray Transfer Screw - Short	24.00	3.10	RP5010
	RP Open Tray Transfer Screw - Long	28.00	3.10	RP5014
	RP Closed Tray Transfer Screw - Short	16.00	2.80	RP5126
	RP Closed Tray Transfer Screw - Long	20.00	2.80	RP5127
	RP Retrieval Screw for Abutments	12.60	2.25	RP3401



Hand & Torque Drivers



Star Hand Driver Short 23.4L RP0064

*For Angeld TMA/TMA prosthetic screw



Star Hand Driver Long 23.4L RP0004

*For Angeld TMA/TMA prosthetic screw



TMA/FC Female Hand Hex Driver - Short 22.4L RS6193

*For straight TMA/Flat connection



TMA/FC Female Hand Hex Driver Long 28.4L RS6191

*For straight TMA/Flat connection



RP 2.2mm Hex Torque Driver Short 20.0L RP0067



RP 2.2mm Hex Torque Driver-Long 25.0L RP0066



Star Torque Driver Short 20.40L RP0065

*For Angeld TMA/TMA prosthetic screw



Star Torque Driver - Long 25.4L RP0005

*For Angeld TMA/TMA prosthetic screw



TMA/FC Hex Torque Driver - Short 20.0L RS6197

*For straight TMA/Flat connection



TMA/FC Hex Torque Driver Long 25.0L RS6196

*For straight TMA/Flat connection



Torque Driver for Ball Attachment 16.7L RS6090



Wrench Ratchet Surgical Instruments Handle 16.70L RA0050



Handpiece Drivers



RP 2.25mm Hex Handpiece Adapter Short 23.7L

RP0068



RP 2.25mm Hex Handpiece Adapter Medium 27.6L

RP0071



RP 2.25mm Hex Handpiece Adapter Long 36.0L

RP0074



Star Handpiece Driver - Short 25.0L

RP0018



Star Handpiece Driver - Long 30.0L

RP0019



RP Handpiece Hex Adapter - Short 20.0L

RP0062

*For Angeld TMA/TMA prosthetic screw

*For Angeld TMA/TMA prosthetic screw



RP Handpiece Hex Adapter - Long 25.0L

RP0072



TMA/FC Hex Handpiece Adapter - Short 20.0L

RS6194

*For straight TMA/Flat connection



TMA/FC Hex Handpiece Adapter - Long 25.0L

RS6195A


*For straight TMA/Flat connection





WP CloseFit™

Wide Platform



Recommended insertion torque up to 50 Ncm

CloseFit™ WP	Ø	L	D3	Catalog No.
	4.3	8.0	3.3	ISPX0843
		10.0	3.2	ISPX1043
		11.5	3.1	ISPX1143
		13.0	3.0	ISPX1343
		15.0	2.8	ISPX1543
		18.0	2.6	ISPX1843
	5	8.0	3.9	ISPX0849
		10.0	3.9	ISPX1049
		11.5	3.8	ISPX1149
		13.0	3.7	ISPX1349
		15.0	3.5	ISPX1549
		18.0	3.3	ISPX1849


Radiographic Template RP CloseFit™ | ASD0612

Drill Sequence Recommended drilling speed up to 1200 rpm

Implant Diameter

ø4.3

ø5.0



You may use Tri-Step drill instead of ø2, ø2.8 and ø3.2 drills in sequence.

Cover Screw & Healing Abutments Recommended Torque 15(Ncm)

WP CloseFit™ Conical Cover Screw




WP0073

WP CloseFit™ Slim Healing Abutment



WP0046, GH2/L7.3
 WP0047, GH3/L8.3
 WP0048, GH4/L9.3
 WP0049, GH5/L10.3
 WP0050, GH6/L11.3

WP CloseFit™ Healing Abutment



WP0022, GH2/L7.3
 WP0023, GH3/L8.3
 WP0024, GH4/L9.3
 WP0025, GH5/L10.3
 WP0026, GH6/L11.3

WP CloseFit™ Conical Healing Abutment




WP0090, GH2/L7.9
 WP0091, GH3/L8.9
 WP0092, GH4/L9.9
 WP0093, GH5/L10.9
 WP0094, GH6/L11.9




WP™ Impression

Recommended Torque 15(Ncm)





Our components for taking impressions for all restoration types can be successfully utilized for digital or conventional protocols either implant-level or abutment level procedures.

Our broad range of solutions and flexible designs enable easy and accurate restorative planning, whether you're dealing with a single tooth, multiple teeth, bridge or full mouth rehabilitation.




Transfers		L	∅	D1	H	GH	RH	Catalog No.
	WP Slim Closed Tray Transfer - Short	11.30	4.00	3.35	9.00	-	-	WP5109
	WP Slim Closed Tray Transfer - Long	15.30	4.00	3.35	13.00	-	-	WP5113
	WP Closed Tray Transfer - Short	11.30	5.20	3.35	9.00	-	-	WP5209
	WP Closed Tray Transfer - Long	15.30	5.20	3.35	13.00	-	-	WP5213
	WP Wide Closed Tray Transfer - Short	11.30	6.00	3.35	9.00	-	-	WP5309
	WP Wide Closed Tray Transfer - Long	15.30	6.00	3.35	13.00	-	-	WP5313
	WP Slim Open Tray Transfer - Short	14.30	6.00	3.35	10.00	-	-	WP5610
	WP Slim Open Tray Transfer - Long	18.30	6.00	3.35	14.00	-	-	WP5614
	WP Open Tray Transfer - Short	14.30	5.20	3.35	10.00	-	-	WP5510
	WP Open Tray Transfer - Long	18.30	5.20	3.35	14.00	-	-	WP5514
	WP Wide Open Tray Transfer - Short	14.30	6.00	3.35	10.00	-	-	WP5410
	WP Wide Open Tray Transfer - Long	18.30	6.00	3.35	14.00	-	-	WP5414
	WP Scan Body	9.3	-	3.35	7.0	-	-	WP0088
		12.3	-	3.35	10.0	-	-	WP0089

Analog / Digital Analog		L	∅	D1	H	GH	RH	Catalog No.
	WP Implant Analog	13.10	4.00	3.35	-	-	-	WP0007
	WP Digital Model Analog 15 (Ncm)	10.00	4.00	3.35	-	-	-	WP1000



Cement Retained Abutments

Recommended Torque 35(Ncm)





Abutments	L	∅	D1	D2	H	GH	RH	Catalog No.
 WP Straight Slim Abutment No End Line	11.00	4.20	3.35	3.36	8.70	2.70	6.00	WP0020
 WP Straight Abutment	9.80	5.05	3.35	3.57	7.50	1.00	6.50	WP0012
	10.80	5.05	3.35	3.57	8.50	2.00	6.50	WP0011
	11.80	5.05	3.35	3.57	9.50	3.00	6.50	WP0010
 WP Angled Abutment 15°	9.80	5.05	3.35	2.61	7.50	1.00	6.50	WP0013
	10.80	5.05	3.35	2.61	8.50	2.00	6.50	WP0014
	11.80	5.05	3.35	2.61	9.50	3.00	6.50	WP0015
 WP Engaging Temporary Abutment	12.65	4.20	3.35	3.90	10.35	0.90	9.45	WP0043
 WP Non-Engaging Temporary Abutment	12.65	4.20	3.35	3.90	10.35	0.90	9.45	WP0044
 WP Plastic Cylinder Abutment (Casting)	12.65	4.20	3.35	3.40	10.35	0.90	9.45	WP0045
 WP Engaging Titanium/Plastic Abutment (Casting)	13.80	5.20	3.35	5.20	11.50	1.50	10.00	WP3095
	13.80	5.20	3.35	5.20	11.50	1.50	10.00	WP3097



TMA™ - Multi Unit

Recommended Torque 35(Ncm)








Our TMA™ (Trans Mucosal Abutment) system is indicated for multiple-unit, screw-retained restorations and may be used in combination with an implant level framework design.

The TMA™ system is used to elevate the restorative seating platform when restoration at implant level is not indicated or practical due to the depth or angle of the implant.

The TMA™ System is designed to allow better prosthetics access when using implants with divergent angles accommodate full and partial edentulous arches, especially when using angled implants and all-on-x technique for full arch restoration.

The TMA™ abutment system is available in straight or angled (17°,30°,45°), engaging the Conical Connection for indexing, with a selection of collar heights. Abutments are delivered with handles for easy insertion and to assist with seating.

TMA		L	Ø	D1	H	GH	RH	Catalog No.
	WP Straight Trans Mucosal Abutment	4.30	4.90	3.35	2.0	1.0	2.10	WP3731
		5.30	4.90	3.35	3.0	2.0	2.10	WP3732
		6.30	4.90	3.35	4.0	3.0	2.10	WP3733
		7.30	4.90	3.35	5.0	4.0	2.10	WP3729
		8.30	4.90	3.35	6.0	5.0	2.10	WP3730
	WP Angled Trans Mucosal Abutment 17° 2.5mm	6.7	4.90	3.35	-	2.54	2.15	WP3734
	WP Angled Trans Mucosal Abutment 17° 3.5mm	7.44	4.90	3.35	-	3.56	2.15	WP3737
	Angled Trans Mucosal Abutment 30° 3mm	6.8	4.90	3.35	-	3.45	2.15	WP3735
	WP Angled Trans Mucosal Abutment 30° 4mm	8.04	4.90	3.35	-	4.11	2.15	WP3738


TMA components		L	∅	D1	D2	H	GH	RH	Catalog No.
	TMA Healing Cap	5.5	4.9	-	3.5	-	-	-	RS5005
		7	4.9	-	3.5	-	-	-	RS5003
	TMA Impression Coping Open Tray	11.0	4.9	-	3.1	-	-	-	RS5006
	TMA Closed Tray Transfer	8.0	4.9	-	4.0	-	-	-	RS5011
	TMA Scan Body	6	4.9	6.5	3.8	4.9	-	-	TM0006
		9	4.9	6.5	3.8	4.9	-	-	TM0009
	TMA Titanium Temporary	12.0	4.9	-	3.5	-	1.5	10.5	RS4900
	TMA Cementing Cone	4.9	4.9	-	2.9	-	1.5	3.4	TM0010
	TMA Plastic Casting Sleeve w/ Prosthetic Screw	12.0	4.9	-	3.3	-	1.5	10.5	RS5001
	TMA Abutment Replica	16.8	4.9	-	1.7	-	-	-	RS5004
	TMA Digital Model Replica Recommended torque 15(Ncm)	12	4.9	-	3.8	-	-	-	RS1001










Overdenture Restoration

Recommended Torque 35(Ncm)



Ball Attachment		L	Ø	D1	D2	H	GH	RH	Catalog No.
	WP Ball Attachment	9.32	3.50	3.35	2.50	4.00	0.50	2.72	WP0051
		9.82	3.50	3.35	2.50	4.50	1.00	2.72	WP0052
		10.82	3.50	3.35	2.50	5.50	2.00	2.72	WP0053
		11.82	3.50	3.35	2.50	6.50	3.00	2.72	WP0054
		12.82	3.50	3.35	2.50	7.50	4.00	2.72	WP0055
		13.82	3.50	3.35	2.50	8.50	5.00	2.72	WP0056



Ball Attachment caps		L	Ø	D1	D2	H	GH	RH	Catalog No.
	Plastic Ball Cap White (Nylon)	2.60	4.10	-	3.60	-	-	-	RS2660
	Plastic Ball Cap Pink Soft (Nylon)	2.60	4.10	-	3.60	-	-	-	RS2662
	Titanium Ball Cap	2.90	4.70	-	4.70	-	-	-	RS2670
	Stainless Steel Ball Cap	2.90	4.70	-	4.70	-	-	-	RS2675SS

Ball Attachment		L	Ø	D1	D2	H	GH	RH	Catalog No.
	Ball Attachment Abutment Replica	16.60	4.00	-	2.50	-	0.80	2.70	RS5740
















Screw - Retained

Recommended Torque 35(Ncm)



Flat Connection	L	Ø	D1	D2	H	GH	RH	Catalog No.
	7.40	4.50	3.35	-	2.0	2.0	-	WP0030
WP Flat Connection Abutment	9.40	4.50	3.35	-	4.0	4.0	-	WP0031
	11.50	4.50	-	3.90	11.10	0.40	10.70	FC5515
	13.60	4.60	-	3.00	13.60	0.10	13.50	FC6011
	5.00	2.50	-	-	-	-	-	FC6015
	3.80	4.50	-	3.10	3.40	0.60	2.80	FC6013
	17.00	-	-	2.50	-	-	-	FC6017
	3.80	4.50	-	3.20	3.40	0.60	2.80	FC6018
	11.00	4.50	-	4.20	-	-	10.2	FC5512
	12.60	4.50	-	3.4	-	-	-	FC5737

Screws		L	Ø	Catalog No.
	WP Abutment Screw	7.5	2.55	WP0001
	TMA Prosthetic Screw 15(Ncm)	4.1	2.1	RS3404
	WP Retaining Screw for TMA	6.50	3.10	WP3406
	WP Retaining Screw for TMA	7.50	3.10	WP3407
	WP Retaining Screw for TMA	8.50	3.10	WP3408
	WP Retaining Screw for TMA	9.50	3.10	WP3409
	WP Retaining Screw for TMA	10.50	3.10	WP3412
	WP Retaining Screw for Angled TMA	7.10	2.35	WP3403
	WP Open Tray Transfer Screw - Short	24.00	3.10	WP5010
	WP Open Tray Transfer Screw - LONG	28.00	3.10	WP5014
	WP Closed Tray Transfer Screw - Short	16.50	3.10	WP5126
	WP Closed Tray Transfer Screw - Long	20.00	3.10	WP5127
	WP Retrieval Screw for Abutments	12.30	2.55	WP3401



Hand & Torque Drivers



Star Hand Driver
Short 23.4L

RP0064

*For Angeld TMA/TMA prosthetic screw



Star Hand Driver
Long 23.4L

RP0004

*For Angeld TMA/TMA prosthetic screw



TMA/FC Female
Hand Hex Driver - Short
22.4L

RS6193

*For straight TMA/Flat connection



TMA/FC Female
Hand Hex Driver
Long 28.4L

RS6191

*For straight TMA/Flat connection



WP 2.6mm
Torque Driver
Short 20.0L

WP0067



WP 2.6mm Torque
Driver-Long 25.0L

WP0066



Star Torque
Driver
Short 20.40L

RP0065

*For Angeld TMA/TMA prosthetic screw



Star Torque
Driver - Long
25.4L

RP0005

*For Angeld TMA/TMA prosthetic screw



TMA/FC Hex Torque
Driver - Short 20.0L

RS6197

*For straight TMA/Flat connection



TMA/FC Hex
Torque Driver
Long 25.0L

RS6196

*For straight TMA/Flat connection



Torque Driver for
Ball Attachment
16.7L

RS6090



Wrench Ratchet
Surgical Instruments
Handle 16.70L

RA0050



Handpiece Drivers



WP 2.65mm
Handpiece Adapter
Short 23.85L

WP0068



WP 2.65mm
Handpiece Adapter
Medium 27.85L

WP0071



WP 2.65mm
Handpiece Adapter
Long 35.95L

WP0074



Star Handpiece
Driver - Short 25.0L

RP0018



Star Handpiece
Driver - Long 30.0L

RP0019



WP Handpiece Hex
Adapter - Short 20.0L

WP0062

*For Angeld TMA/TMA prosthetic screw

*For Angeld TMA/TMA prosthetic screw



WP Handpiece Hex
Adapter - Long 25.0L

WP0072



TMA/FC Hex
Handpiece
Adapter - Short 20.0L

RS6194

*For straight TMA/Flat connection



TMA/FC Hex Handpiece
Adapter - Long 25.0L

RS6195A

*For straight TMA/Flat connection







Tri-Step & Twist Step Surgical Drills



Our DLC (Diamond Like Carbon) Tri-Step/Twist Step Surgical Drills are uniquely designed and coated to improve the surgical drilling process.

Tri-Step & Twist-Step DLC drills simplify the cutting process, save time and minimize risks.

Our improved DLC drill coating provides state-of-the-art hardness and slickness that increases lifespan of the drill, maintains sharpness for extended drilling and generates less heat, thereby reducing chances of heating the bone or complications arising from conventional drilling protocols.

The drills come in the following diameters: 2.0, 2.5, 2.8, 3.2, 3.6, 4.2, 4.6, 5.2, 5.6

Specially designed for Adin's implants and suitable for their simple and convenient drilling protocol. They are clearly and conveniently marked and can be used with or without drill stops.



Tri-Step drills





















Adin's Tri-Step drills condense drilling protocols and avoid the need for using multiple drills throughout complex procedures.

Because friction at the cutting site is reduced, complications common with high speed drilling are decreased and faster & better recovery can occur.

L	Ø		Size	Catalog No.
	1.8-3.2		short	TD1511
	1.8-3.2		long	TD2218








Twist Step drills

L	Ø		Color	Catalog No.
short  6-13mm	2.0		white	DL1320
	1.8-2.5		purple	DL1325
	2.3-2.8		red	DL1328
	2.8-3.2		blue	DL1332
	3.2-3.6		green	DL1336
	3.8-4.2		black	DL1342
	4.2-4.6		yellow	DL1346
	4.8-5.2		brown	DL1352
	5.2-5.6		grey	DL1356
	long  8-18mm	2.0		white
1.8-2.5			purple	DL1825
2.3-2.8			red	DL1828
2.8-3.2			blue	DL1832
3.2-3.6			green	DL1836
3.8-4.2			black	DL1842
4.2-4.6			yellow	DL1846
4.8-5.2			brown	DL1852
5.2-5.6			grey	DL1856



Tools



		Catalog No.
	Surgical Torque Ratchet 35-100 Ncm	RS6111
	Prosthetic Torque Ratchet 15-35 Ncm	RS6112
	Wrench Ratchet	RS6100
	Surgical Driver	RS7001



Depth Probe

RS6180

Tools



Parallel Pin 11mm

RS6150



Parallel Pin 16mm

RS6155



Drill Extension


RS9027

Osteotomes



Osteotome Handle

HA0074



Osteotome Insert 2.0mmD

HA0061

Osteotome Insert 2.8mmD

HA0064

Osteotome Insert 3.25mmD

HA0067

Osteotome Insert 3.65mmD


HA0070







Surgical Kits Conical connection



Successful surgical procedures start with having the right instruments and tools to ensure your desired treatment outcomes are met.

That's why our range has everything you need, including drills, keys, ratchets and insertion tools. Not only do they make procedures safer and simpler, they result in an enhanced patient experience as well.

Our Surgical Kits streamline and simplify your implant procedures.

Surgical Kits contain			Premium CF6300	Advanced CF6302	Starter CF6303
Drills:					
	RS8032	Tip Drill	V		
	TD2218	Tri-Step Drill-Long	V		
	DL1820	Pilot Drill 2mmD x 8-18mmL	V	V	V
	DL1825	Twist Step Drill 1.8/2.5mmD x 8-18mmL	*		
	DL1828	Twist Step Drill 2.3/2.8mmD x 8-18mmL	V	V	V
	DL1832	Twist Step Drill 2.8/3.2mmD x 8-18mmL	V	V	V
	DL1836	Twist Step Drill 3.2/3.6mmD x 8-18mmL	V	V	V
	DL1842	Twist Step Drill 3.8/4.2mmD x 8-18mmL	V	V	
	DL1846	Twist Step Drill 4.2/4.6mmD x 8-18mmL	V	V	
Drills:					
	TD1511	Tri-Step Drill - Short	V		
	DL1320	Pilot Drill 2mmD x 6-13mmL	V		
	DL1325	Twist Step Drill 1.8/2.5mmD x 6-13mmL	*		
	DL1328	Twist Step Drill 2.3/2.8mmDx6-13mmL	V		
	DL1332	Twist Step Drill 2.8/3.2mmDx6-13mmL	V		
	DL1336	Twist Step Drill 3.2/3.6mmDx6-13mmL	V		
	DL1342	Twist Step Drill 3.8/4.2mmDx6-13mmL	V		
	DL1346	Twist Step Drill 4.2/4.6mmDx6-13mmL	V		
Surgical tools:					
	RS6111	Surgical Torque Ratchet 35-100 Ncm	V	V	V
	RA0050	Wrench Ratchet Surgical Instruments Handle	V	V	V
	RS6180	Depth Probe	V	V	
Parallel pins:					
	RS6150	Parallel Pin Short (qty. 4)	V	V	
	RS6155	Parallel Pin Long (qty. 2)	V	V	
Insertion tools:					
	NP0068	NP 2mm Hex Handpiece Adapter - Short	V	V	V
	NP0071	NP 2mm Hex Handpiece Adapter - Medium	V	V	V
	NP0074	NP 2mm Hex Handpiece Adapter - Long	V	V	V
	RP0068	RP 2.25mm Hex Handpiece Adapter - Short	V	V	V
	RP0071	RP 2.25mm Hex Handpiece Adapter - Medium	V	V	V
	RP0074	RP 2.25mm Hex Handpiece Adapter - Long	V	V	V
	WP0068	WP 2.65mm Hex Handpiece Adapter - Short	V	V	V
	WP0071	WP 2.65mm Hex Handpiece Adapter - Medium	V	V	V
	WP0074	WP 2.65mm Hex Handpiece Adapter - Long	V	V	V
Drivers:					
	RP0064	Star Hand Driver - Short	V	V	V
	RP0004	Star Hand Driver - Long	V	V	V
Drill extension:					
	RS9027	Drill Extension	V	V	

*Products may be ordered separately as an addition to the kit

Premium Kit CF6300



Drills	Catalog No.
Tip Drill	RS8032
Tri-Step Drill-Long	TD2218
Pilot Drill 2mmD x 8-18mmL	DL1820
Twist Step Drill 1.8/2.5mmD x 8-18mmL	DL1825
Twist Step Drill 2.3/2.8mmD x 8-18mmL	DL1828
Twist Step Drill 2.8/3.2mmD x 8-18mmL	DL1832
Twist Step Drill 3.2/3.6mmD x 8-18mmL	DL1836
Twist Step Drill 3.8/4.2mmD x 8-18mmL	DL1842
Twist Step Drill 4.2/4.6mmD x 8-18mmL	DL1846
Tri-Step Drill - Short	TD1511
Pilot Drill 2mmD x 6-13mmL	DL1320
*Twist Step Drill 1.8/2.5mmDx6-13mmL	DL1325
Twist Step Drill 2.3/2.8mmDx6-13mmL	DL1328
Twist Step Drill 2.8/3.2mmDx6-13mmL	DL1332
Twist Step Drill 3.2/3.6mmDx6-13mmL	DL1336
Twist Step Drill 3.8/4.2mmDx6-13mmL	DL1342
Twist Step Drill 4.2/4.6mmDx6-13mmL	DL1346

Surgical Tools	Catalog No.
Wrench Ratchet Surgical Instruments Handle	RA0050
Surgical Torque Ratchet 35-100 Ncm	RS6111
Depth Probe	RS6180

Paraller Pin	Catalog No.
Parallel Pin Short 11mm (qty. 4)	RS6150
Parallel Pin Long 16mm (qty. 2)	RS6155

Insersion Tools	Catalog No.
NP 2mm Hex Handpiece Adapter - Short	NP0068
NP 2mm Hex Handpiece Adapter - Medium	NP0071
NP 2mm Hex Handpiece Adapter - Long	NP0074
RP 2.25mm Hex Handpiece Adapter - Short	RP0068
RP 2.25mm Hex Handpiece Adapter - Medium	RP0071
RP 2.25mm Hex Handpiece Adapter - Long	RP0074
WP 2.65mm Hex Handpiece Adapter - Short	WP0068
WP 2.65mm Hex Handpiece Adapter - Medium	WP0071
WP 2.65mm Hex Handpiece Adapter - Long	WP0074

Drivers	Catalog No.
Star Hand Driver - Short	RP0064
Star Hand Driver - Long	RP0004

Drill Extension	Catalog No.
Drill Extension	RS9027

*Products may be ordered separately as an addition to the kit

Advanced Kit CF6302



Drills	Catalog No.
Pilot Drill 2mmD x 8-18mmL	DL1820
Twist Step Drill 2.3/2.8mmD x 8-18mmL	DL1828
Twist Step Drill 2.8/3.2mmD x 8-18mmL	DL1832
Twist Step Drill 3.2/3.6mmD x 8-18mmL	DL1836
Twist Step Drill 3.6/4.2mmD x 8-18mmL	DL1842
Twist Step Drill 4.2/4.6mmD x 8-18mmL	DL1846

Insertion Tools	Catalog No.
NP 2mm Hex Handpiece Adapter - Short	NP0068
NP 2mm Hex Handpiece Adapter - Medium	NP0071
NP 2mm Hex Handpiece Adapter - Long	NP0074
RP 2.25mm Hex Handpiece Adapter - Short	RP0068
RP 2.25mm Hex Handpiece Adapter - Medium	RP0071
RP 2.25mm Hex Handpiece Adapter - Long	RP0074
WP 2.65mm Hex Handpiece Adapter - Short	WP0068
WP 2.65mm Hex Handpiece Adapter - Medium	WP0071
WP 2.65mm Hex Handpiece Adapter - Long	WP0074

Surgical Tools	Catalog No.
Wrench Ratchet Surgical Instruments Handle	RA0050
Surgical Torque Ratchet 35-100 Ncm	RS6111
Depth Probe	RS6180


Drivers	Catalog No.
Star Hand Driver - Short	RP0064
Star Hand Driver - Long	RP0004

Parallel Pin	Catalog No.
Parallel Pin Short 11mm (qty. 4)	RS6150
Parallel Pin Long 16mm (qty. 2)	RS6155

Drill Extension	Catalog No.
Drill Extension	RS9027

*The advanced kit is also available with short drills | Catalog No. CF6304

Starter Kit CF6303



Drills	Catalog No.
Pilot Drill 2mmD x 8-18mmL	DL1820
Twist Step Drill 2.3/2.8mmD x 8-18mmL	DL1828
Twist Step Drill 2.8/3.2mmD x 8-18mmL	DL1832
Twist Step Drill 3.2/3.6mmD x 8-18mmL	DL1836
Twist Step Drill 3.6/4.2mmD x 8-18mmL	DL1842

Insertion Tools	Catalog No.
NP 2mm Hex Handpiece Adapter - Short	NP0068
NP 2mm Hex Handpiece Adapter - Medium	NP0071
NP 2mm Hex Handpiece Adapter - Long	NP0074
RP 2.25mm Hex Handpiece Adapter - Short	RP0068
RP 2.25mm Hex Handpiece Adapter - Medium	RP0071
RP 2.25mm Hex Handpiece Adapter - Long	RP0074
WP 2.65mm Hex Handpiece Adapter - Short	WP0068
WP 2.65mm Hex Handpiece Adapter - Medium	WP0071
WP 2.65mm Hex Handpiece Adapter - Long	WP0074

Surgical Tools	Catalog No.
Wrench Ratchet Surgical Instruments Handle	RA0050
Wrench Ratchet	RS6100

Drivers	Catalog No.
Star Hand Driver - Short	RP0064
Star Hand Driver - Long	RP0004



Founded in 2001, Adin Dental Implant Systems is a global leader in the design, development and manufacture of dental implant solutions & on track to become the world's leading independent dental implant provider.

Every business starts with a vision. Ours is simple:
We want to create the best possible dental implant solutions
that offer uncompromising quality at an affordable price.
It's our people that help our business thrive and grow.
That's why we place so much importance on building strong
personal relationships that enable us to meet our doctors'
and distributors' needs. Understanding that their success is
our success, we're focused on providing professional, high
quality and affordable solutions and exceptional service
that help them grow.

