

Take your practice to the next level











# | Guided Surgery Approach

Using Guided surgery for implant placement has become increasingly popular among dentists.

There are many reasons to use guided surgery:


#### **Patient Benefits**

- Surgical time reduction ideal for patient comfort and apprehension
- Easy to execute flapless surgery significantly improves healing process and time in majority of cases
- Minimization of surgical complications

#### **Dentist Benefits:**

- Significantly improves placement accuracy
- Significantly reduces likelihood of complications
- Minimizes patient chair time efficient utilization of time


# | Why choose adinguide™?




The adinguide™ was planned with you in mind, bringing the guided surgery world to your clinic. Specifically designed to support effortless procedures, accurate and predictable implant placement.

# | Keyless Procedure

This innovative approach developed to support effortless procedures, accurate and predictable implant placement.




Other solutions .....



## | ActiveFlow™ Irrigation Technology


The adinguide<sup>™</sup> ActiveFlow<sup>™</sup> Irrigation Technology is a unique design that forces coolant through the guide and ensure that irrigation reaches the bone. Within each cylinder are unique channels which allow coolant to stream through the sleeve while drilling.

- $\bullet$  ActiveFlow  $\protect\ensuremath{^{\text{T}}}$  Technology through the guide irrigation
- As the drill starts to spin, ActiveFlow™ starts operation
- Delivers the coolant from the hand piece directly to the surgical site and may prevent bone heating



# | Self-Centering Drills & Tools

Our keyless, self - centering drills and tools are designed to enable dentists to perform faster, easier and more intuitive implant surgeries. Organized in a way that enhances safety and ease of use, Adin sleeves are embedded within the 3D printed surgical guide and precisely center and guide drills and insertion tools. Self - centering drills, controlled by a main diameter cylinder at the top of each drill, centers the drill in the guide's sleeve. The surgical guide and drilling procedure are designed so that the self-centering cylinder engages the sleeve before the drill touches the bone in order to ensure correct positioning of the drill.



### | Built-in Stoppers Enhance Safety

Built-in stoppers ensure precise and accurate drilling to the desired depth.

The fixed value of the offset between the top of the implant and sleeves allow doctors to decide to use shorter drills than planned. Enabling another safely measure, allowing another layer of verification before reaching the final drilling depth.



# | Simplicity

The adinguide™ kits are Intuitive, simple and effortless:

No need for keys, or any other assisting tools. Hassle free, logical and straightforward procedures.



# Guided Surgery Kit RS Touareg™, CloseFit™ RP/WP Implants Drills GSK001



# | Indicators for Ease of Use

Procedure advances from left to right, features Intelligent and logical drill layout

- ullet Each column is marked and color coded for a specific drill and implant diameter  $oldsymbol{arnothing}$
- Color coding with simple direction
- · Each row is marked for a specific length
- Indication arrows for direction of workflow

#### | Drill Indications

Each drill is color coded according to drill diameter and also laser marked for the drill length and diameter. The diameter of each drill group refers to the diameter of the center of the drill, as indicated in the drawing on the right side of the tray. The drawing on the left side provides easy verification of the length of the drills by placing the large cylinders in the carved out area.





#### | Tissue Punch

The tissue punch is utilized for cutting through soft tissue. It is designed to make a circular incision of 5mm around each planned implant site. The tissue punch is designed to self-center to the surgical guide. It is recommended to use 30 rpm. After punching, take the guide off and remove soft tissue with a compatible instrument.



#### | Pilot Drill

The starter pilot drill is used to make the initial penetration into bone. It is design for optimal initial centering and also has a reamer that can flatten bone, where applicable.



### | Spare Space

This space can be utilized for 18mm drills or any other use.


## | Implant Insertion


Adin's Guided Surgery System is designed to allow insertion of mountless implants through the sleeves of the surgical guide.

This ensures accurate 3D positioning of the implants.

The insertion tools allow rigid connection to mountless implants by means of a holding screw, delivery and insertion of the implant to its correct position, using either the hand operated ratchet or the handpiece adaptor.

Recommended insertion torque up to 50 Ncm. To ensure correct orientation of the Hex, the final position of the implant may be adjusted after removal of the surgical guide.









# O | Guided Surgery Kit





#### For full guided procedure:





RS Touareg™ S/OS Implant Platform



GS Tools Kit RS GSK003











GS Accessories KIT - RP&WP GSK004







GS1025 Drill S-5.5, 2.5/1.8mmDx10.0mmL



GS1125 Drill S-5.5, 2.5/1.8mmDx11.5mmL



GS1325 Drill S-5.5, 2.5/1.8mmDx13.0mmL



GS1625 Drill S-5.5, 2.5/1.8mmDx16.0mmL





16 Ø3.2 L16

GS1632 Drill S-5.5, 3.2/2.5mmDx16.0mmL

**G\$5055** Tissue Punch 5.0mmD, S-5.5

GS0625 Pilot Drill S-5.5, 2.5/1.8mmDx6.0mmL

Ø2.5 L6













GS1346 Drill S-5.5, 4.6/3.85mmDx13.0mmL



GS1646 Drill S-5.5, 4.6/3.85mmDx16.0mmL



# Guided Surgery Kit CloseFit™ UNP/NP Narrow Implants Drills & Tools GSK002




For full guided procedure:














RP0004

Star Hand Driver - Long

16 GS0023

Manual Socket L10.5

# © Guided Surgery Kit RS Touareg™ Implants Tools GSK003



For full guided procedure:















RS Grip Fixture Mount

**RS Extraction Tool** 





GS1915

Anchoring Drill S-2.0, 2.0/1.5mmDx19mmL

Ø1.5 L19



SPARE

SURGICAL DRIVER

TORQUE RATCHET

TORQUE RATCHET

TORQUE RATCHET

TORQUE RATCHET

TORQUE RATCHET

EXTRACTION

EXTRACTION

Gental implant solutions

Guided Surgery Tools Kit

RS Implants

DRIVERS

EXTENTION

HEX DRIVER

TORQUE RATCHET

TORQUE RATCHE



Manual Socket L14.5

# © | Guided Surgery Kit CloseFit™ RP/WP Implants Tools GSK004



For full guided procedure:







GS Drill Kit GSK001







Anchoring Drill S-2.0, 2.0/1.7mmDx19mmL



Handpiece Guide Fixation Screw Driver



Guided RP Ratchet Fixture Mount S-5.5, Screw Retained



RP Ratchet retrieval tool



Guided RP Handpiece Fixture Mount S-5.5



Guided RP Ratchet Fixture Mount S-5.5




Guided WP Handpiece Fixture Mount S-5.5




Guided WP Ratchet Fixture Mount S-5.5



Guided WP Ratchet Fixture Mount S-5.5, Screw Retained







#### | Purchase Options

Catalog No.	Description	Package contents	
GS0006	Guided Vertical Fixation Screw S-3.3		x 3 pack
GS0007	Guided Vertical Fixation Screw S-5.5		x 3 pack
G\$1920	2.0mmD Anchoring Screw		x 3 pack
GS0055	5.5mmD Sleeve		x 3 pack
			x 10 pack
GS0033	3.3mmD Narrow Implant Sleeve		x 3 pack
			x 10 pack
GS0020	2.0mmD Anchoring Sleeve		x 3 pack



#### 18mm Drills



Catalog No. GS1825



Drill S-5.5, 3.2/2.5mmDx18.0mmL

Catalog No. GS1832



Catalog No. GS1838



Drill S-5.5, 4.6/3.85mmDx18.0mmL

Catalog No. GS1846



Drill S-3.3, 2.5/1.8mmDx18.0mmL


Catalog No. GSN1825




Illustrations are not to scale

Availability of the products in some places is limited and depends on local regulations.

For more information, please contact your local representative.









#### **About Adin**

Every business starts with a vision.

Ours is simple: We want to create the best possible dental implant solutions that offer uncompromising quality at an affordable price. It's our people that help our business thrive and grow. That's why we place so much importance on building strong personal relationships that enable us to meet our doctors and distributors' needs. Understanding that their success is our success, we're focused on providing professional, high quality and affordable solutions and exceptional service that help them grow.



