




**A New Standard in
Endodontic Treatment**

Effective · Simple · Gentle



Distributed by:

Adin Dental Implant Systems Ltd.
Industrial Zone Alon Tavor,
POB 1128, Afula 1811101, Israel
gentlefile@adin-implants.com
www.adin-implants.com

Manufactured by:

MedicNRG Ltd.
Kibbutz Afikim 15148
Israel
www.gentlefile.com



Effective · Simple · Gentle

With an innovative filing technique and unique patented design geared to protect the canal's anatomical structure and avoid dentine perforation, Gentlefile helps every dentist perform optimal root canal treatments.

Optimal Root Canal Treatment

Continuous flexibility throughout file length

- Adapts to any canal shape and maintains its anatomical structure
- Prevents dentine perforations, cracks and ledges
- Optimal access to posterior teeth

Automatic rotational speed with centrifugal effect


- Quickly and gently files the canal's full surface area
- Follows the canal's cross sectional shape in its entirety
- Activates the antibacterial fluid along the canal's entire surface and length
- Significantly removes the smear layer and improves access to lateral canals
- Cleans canals with irregular cross-sections




Continuous flexibility throughout file length



Access to irregularities in the canal




Surface after treatment-
the Coronal Third
Enlargement: X 2500



A combination of file flexibility and centrifugal movement permits accessing and cleaning of the irregular parts of canals that NiTi files cannot reach.

Safe Root Canal Treatment

- Files canal walls with no drilling or cutting
- Prevents torque overload
- Minimizes file breakage rate



When compared to NiTi files, Gentlefile[®] generates far less force on the canal wall, therefore causing no damage to the dentin.

Simple Root Canal Treatment

- Simple, quick preparation process with minimal files quantity (2-3)
- One file length for all canal shapes
- Each file is designed with a different ratio between vertical force and flexibility, adapted for all treatments
- No need for adjustment, easy to control

