




Adin's long implant range serves as an extension to the current Touareg[™]-OS and Touareg CloseFit[™] WP platform implants, allowing implant placement options up to 25mm in length.

In cases where more bone-to-implant contact is required, Long implants may serve as a solution where the space allows for longer implants, alternatively, long implants can be used effectively as bi-cortical anchorage to compensate for lack of bone while still achieving implant stability in cases where implants of regular length might not.

Adin's long implant range carry specific features designed to enhance their performance such as narrow apex for limited spaces and a smooth neck to better suite gingival engagement.

Adin's unique and reliable Osseofix[™] surface treatment is utilized for both current Touareg[™]-OS and CloseFit[™] WP platform long implants Compatible drill set and depth probe have been designed to maximize ease of use and safety.




Ømm	Lmm	D3 mm	Catalog No.
4.25	20.00	2.50	ISF2042
	22.50	2.50	ISF2242
	25.00	2.50	ISF2542

Touareg[™]-OS Recommended insertion torque up to 50 Ncm


CloseFit[™] WP Recommended insertion torque up to 50 Ncm

Ømm	Lmm	D3 mm	Catalog No.
4.30	20.00	2.50	ISPX2043
	22.50	2.50	ISPX2243
	25.00	2.50	ISPX2543

Extra-long Drills




Tool




Drill Sequence

Recommended drilling speed up to 1200 rpm


Touareg[™]-OS





CloseFit[™] WP







About Adin

Every business starts with a vision.

Ours is simple: We want to create the best possible dental implant solutions that offer uncompromising quality at an affordable price. It's our people that help our business thrive and grow. That's why we place so much importance on building strong personal relationships that enable us to meet our doctors and distributors' needs. Understanding that their success is our success, we're focused on providing professional, high quality and affordable solutions and exceptional service that help them grow.



MB049-1 03.2022



Industrial Park Allon Tavor, POB 1128, Afula 1811101, Israel T. +972-4-6426732 | info@adin-implants.com | www.adin-implants.com