



CloseFit™

Regular and
Wide Implants

Product Catalog



CloseFit™

Regular and Wide Implants

Our strong CloseFit™ Conical Connection is designed for optimal seating while reducing and virtually eliminating micromovements.

Appropriate for immediate function/loading capability in all bone types.

CloseFit™ implants are designed to be used in all indications and in all regions of the lower and upper jaws.

The back-tapered coronal portion is designed to reduce pressure and enable healthy cortical bone around the implant collar

Implant positioning, crestal and subcrestal, are well-documented placement protocols supported by clinical evidence placement protocols.

The use of either option is decided upon by the clinician, and should be tailored to patient-specific factors such as bone size and quality, soft tissue characteristics, functional and esthetics. Adin bone level implants support both methods. Please note that when using the AdinGuide™ system, subcrestal placement requires drilling adaptations, and freehand final positioning after guide removal.

CloseFit™

Regular and Wide

RP CloseFit™

Ø 3.5mmD
Regular Platform

*Built-In Platform Switching

*OsseoFix™ surface




WP CloseFit™

Ø 4.3/5.0mmD
Wide Platform

*Back tapered collar

***Micro Threads:**
Engineered to enhance bone
regeneration, prevent bone
loss & increase bone implant
contact

*Bone - condensing property
*Self - drilling




*The features related to UNP & NP Implants

RP CloseFit™ Regular Platform


Recommended insertion torque up to 50 Ncm

	D1	L	D2	Catalog No.
	3.5	8.0	1.5	ISPX0835
		10.0	1.5	ISPX1035
		11.5	1.5	ISPX1135
		13.0	1.5	ISPX1335
		15.0	1.5	ISPX1535
		18.0	1.5	ISPX1835



Radiographic Template RP CloseFit™ | ASD0612

Drill Sequence for RP Recommended drilling speed up to 1200 rpm



All measurements are in mm

*Products may be ordered separately as an addition to the kit.


Note: (x.xx) Drill to the depth of the cortex only. All measurements are in mm.


⚠ CAUTION: The drill preparation is up to 1mm longer than the implant.

* The drill apical step = 1.8mmD.

WP CloseFit™ Wide Platform

Recommended insertion torque up to 50 Ncm

	D1	L	D2	Catalog No.
	4.3	8.0	1.9	ISPX0843
		10.0	1.9	ISPX1043
		11.5	1.9	ISPX1143
		13.0	1.9	ISPX1343
		15.0	1.9	ISPX1543
		18.0	1.9	ISPX1843
	5	8.0	2.5	ISPX0849
		10.0	2.5	ISPX1049
		11.5	2.5	ISPX1149
		13.0	2.5	ISPX1349
		15.0	2.5	ISPX1549
		18.0	2.5	ISPX1849



Radiographic Template WP CloseFit™ | ASD0611

Drill Sequence for WP Recommended drilling speed up to 1200 rpm

Implant Diameter	ø4.3		ø5.0	
Hard bone				
	ø2.0	ø2.8	ø3.2	*Tri-step
Soft bone				
	ø2.0	ø2.8	ø3.2	*Tri-step

All measurements are in mm

*Products may be ordered separately as an addition to the kit.

Note: (x.xx) Drill to the depth of the cortex only. All measurements are in mm.

⚠ CAUTION: The drill preparation is up to 1mm longer than the implant.

* The drill apical step = 1.8mmD.

Cover Screw & Healing Abutments

Recommended Torque 15(Ncm)



Cover Screw Packed with implant

		Description	Catalog No.
		RP CloseFit™ Conical Cover Screw	RP0073
		WP CloseFit™ Conical Cover Screw	WP0073

Healing Abutments


		Description	GH	Catalog No.
		RP CloseFit™ Slim Healing Abutment	2	RP0046
			3	RP0047
			4	RP0048
			5	RP0049
			6	RP0050
		RP CloseFit™ Healing Abutment	2	RP0022
			3	RP0023
			4	RP0024
			5	RP0025
			6	RP0026
		RP CloseFit™ Conical Healing Abutment	2	RP0090
			3	RP0091
			4	RP0092
			5	RP0093
			6	RP0094
		Description	GH	Catalog No.
		WP CloseFit™ Slim Healing Abutment	2	WP0046
			3	WP0047
			4	WP0048
			5	WP0049
			6	WP0050
		WP CloseFit™ Healing Abutment	2	WP0022
			3	WP0023
			4	WP0024
			5	WP0025
			6	WP0026
		WP CloseFit™ Conical Healing Abutment	2	WP0090
			3	WP0091
			4	WP0092
			5	WP0093
			6	WP0094


Impression Recommended Torque 15(Ncm)

Our components for taking impressions for all restoration types can be successfully utilized for digital or conventional protocols either implant-level or abutment level procedures.

Our broad range of solutions and flexible designs enable easy and accurate restorative planning, whether you're dealing with a single tooth, multiple teeth, bridge or full mouth rehabilitation.

RP Transfers | Recommended Torque 15(Ncm)

	Open Tray	Description	L	H	Catalog No.
		RP Open Tray Transfer Short	19.5	11.8	RP5613
		RP Open Tray Transfer Long	23.5	14.8	RP5617
		RP Open Tray Transfer Screw Short			RP5010
		RP Open Tray Transfer Screw Long			RP5014

	Close Tray	Description	L	H	Catalog No.
		RP Close Tray Transfer Short	10.7	8.8	RP5111
		RP Close Tray Transfer Long	14.9	13.0	RP5115
		RP Open Tray Transfer Screw Short			RP5126
		RP Open Tray Transfer Screw Long			RP5127


RP Scan Bodies


	Description	H	Catalog No.		Description	H	Catalog No.
	RP Scan Body	7	RP0088		RP Scan Body	10	RP0089

RP Analogs

	Description	Catalog No.		Description	Catalog No.
	RP Digital Analog	RP1000		RP Implant Analog	RP0007
	Retaining Screw for Digital Model Analog	RS1002			

WP Transfers | Recommended Torque 15(Ncm)

	Open Tray	Description	L	H	Catalog No.
		WP Open Tray Transfer Short	19.5	11.8	WP5613
		WP Open Tray Transfer Long	23.5	14.8	WP5617
		WP Open Tray Transfer Screw Short			WP5010
		WP Open Tray Transfer Screw Long			WP5014

	Close Tray	Description	L	H	Catalog No.
		WP Close Tray Transfer Short	10.7	8.8	WP5111
		WP Close Tray Transfer Long	15.0	13.0	WP5115
		WP Open Tray Transfer Screw Short			WP5126
		WP Open Tray Transfer Screw Long			WP5127

WP Scan Bodies

	Description	H	Catalog No.		Description	H	Catalog No.
	WP Scan Body	7	WP0088		WP Scan Body	10	WP0089


WP Analogs


	Description	Catalog No.		Description	Catalog No.
	WP Digital Analog	WP1000		WP Implant Analog	WP0007
	Retaining Screw for Digital Model Analog	RS1002			

Temporary & Cement Retained Restoration

Recommended Torque 25(Ncm)

RP Abutments

		Description	D1	GH	RH	Catalog No.
		RP Engaging Temporary Abutment	3.7	1.0	9.5	RP0043
		RP Non-Engaging Temporary Abutment	3.7	1.0	9.5	RP0044
		RP Straight Slim Abutment No End Line	4.0	2.8	6.0	RP0020
		RP Esthetic Straight Abutment	4.4	1.0	6.5	RP0012
			4.4	2.0	6.5	RP0011
			4.4	3.0	6.5	RP0010
		RP Esthetic Angled Abutment 15°	4.4	1.0	6.5	RP0013
			4.4	2.0	6.5	RP0014
			4.4	3.0	6.5	RP0015
		RP Angled Titanium Abutment 25°	4.5	1.0	6.1	RP0251
			4.5	2.0	6.1	RP0252
			4.5	3.0	6.1	RP0253
		RP Angled Titanium Abutment 35°	4.5	1.0	6.1	RP0351
			4.5	2.0	6.1	RP0352
			4.5	3.0	6.1	RP0353
		RP Engaging Plastic Cylinder Abutment (Casting)	3.7	1.0	9.5	RP0045
		RP Engaging Cobalt Chrome/Plastic Abutment (Casting)	4.7	1.5	10.0	RP3094
		RP Engaging Titanium/Plastic Abutment (Casting)	4.7	1.5	10.0	RP3095
		RP Non-Engaging Titanium/Plastic Abutment (Casting)	4.7	1.4	10.0	RP3097
		RP Non-Engaging Cobalt Chrome/Plastic Abutment (Casting)	4.7	1.4	10.0	RP3100
		RP Abutment Screw				RP0001



WP Abutments

		Description	D1	GH	RH	Catalog No.
		WP Engaging Temporary Abutment	4.2	0.9	9.5	WP0043
		WP Non-Engaging Temporary Abutment	4.2	0.9	9.5	WP0044
		WP Straight Slim Abutment No End Line	4.2	2.7	6.0	WP0020
		WP Straight Abutment	5.0	1.0	6.5	WP0012
			5.0	2.0	6.5	WP0011
			5.0	3.0	6.5	WP0010
		WP Angled Abutment 15°	5.0	1.0	6.5	WP0013
			5.0	2.0	6.5	WP0014
			5.0	3.0	6.5	WP0015
		WP Angled Titanium Abutment 25°	5.0	1.0	6.0	WP0251
			5.0	2.0	6.0	WP0252
			5.0	3.0	6.0	WP0253
		WP Angled Titanium Abutment 35°	5.0	1.0	6.0	WP0351
			5.0	2.0	6.0	WP0352
			5.0	3.0	6.0	WP0353
		WP Plastic Cylinder Abutment (Casting)	4.2	0.9	9.45	WP0045
		WP Engaging Titanium/Plastic Abutment (Casting)	5.2	1.5	10.0	WP3095
		WP Non-Engaging Titanium/ Plastic Abutment (Casting)	5.2	1.5	10.0	WP3097
		WP Abutment Screw				WP0001

IRIS™ Snap-on

RP IRIS™ Kit



Description	RH	GH	Catalog No.
IRIS™ RP CloseFit™ Kit	4	1	IR3041
		2	IR3042
		3	IR3043
		4	IR3044
	6	1	IR3061
		2	IR3062
		3	IR3063
		4	IR3064
	8	1	IR3081
		2	IR3082
		3	IR3083
		4	IR3084

RP Snap-on Impression Kit



Description	RH	GH	Catalog No.
IRIS™ RP Snap-on Impression Kit	4	1	IR3141
		2	IR3142
		3	IR3143
		4	IR3144
	6	2	IR3262
		4	IR3264

RP Snap-on Impression Pack



Description	RH	GH	Catalog No.
IRIS™ RP Abutment	4	1	IR0341
		2	IR0342
		3	IR0343
		4	IR0344
	6	1	IR0361
		2	IR0362
		3	IR0363
		4	IR0364
	8	1	IR0381
		2	IR0382
		3	IR0383
		4	IR0384

WP IRIS™ Kit



Description	RH	GH	Catalog No.
IRIS™ WP CloseFit™ Kit	4	1	IR5041
		2	IR5042
		3	IR5043
		4	IR5044
	6	1	IR5061
		2	IR5062
		3	IR5063
		4	IR5064

WP Snap-on Impression Kit



Description	RH	GH	Catalog No.
IRIS™ WP Snap-on Impression Kit	4	1	IR5141
		2	IR5142
		3	IR5143
		4	IR5144
	6	2	IR5262
		4	IR5264

WP Snap-on Impression Pack



Description	RH	GH	Catalog No.
IRIS™ WP Abutment	4	1	IR0541
		2	IR0542
		3	IR0543
		4	IR0544
	6	1	IR0561
		2	IR0562
		3	IR0563
		4	IR0564

Purchase Options

	Description	RH	Catalog No.
	RP Implant Analog		RP0007
	WP Implant Analog		WP0007
	IRIS™ Abutment replica	4	IR0404
		6	IR0406
		8	IR0408



	Description	Qty	Catalog No.
	IRIS™ PEEK Comfort cap	4	IR0180
		6	IR0181
		8	IR0089
	IRIS™, Burnout anti rotation cap red	1	IR0182
		4	IR0185
	IRIS™, Burnout cap white	4	IR0183
		6	IR0186
	IRIS™, Impression coping	1	IR0179
		4	IR0184


Single TMA™ Recommended Torque 35(Ncm)

Adin's STMA offers the possibility to combine the benefits of trans-gingival approach and screw-retained single tooth restoration yielding greater flexibility in both esthetics, safety and restoration maintenance.





RP Abutments

		Description	GH	Catalog No.
		RP Straight Single TMA Abutment	1.0	RP0032S
			2.0	RP0033S
			3.0	RP0034S
			4.0	RP3729S
			5.0	RP3730S
		RP Retaining Screw for TMA 1mm		RP0039
		RP Retaining Screw for TMA 2mm		RP0040
		RP Retaining Screw for TMA 3mm		RP0041
		RP Retaining Screw for TMA 4mm		RP3409
		RP Retaining Screw for TMA 5mm		RP3412

		Description	GH1	GH2	Catalog No.
		RP Angled Single TMA Abutment 17°	1.2	2.6	RP3734S
			2.1	3.5	RP3737S
		RP Angled Single TMA Abutment 30°	1.2	3.7	RP3735S
			1.6	4.0	RP3738S
		RP Retaining Screw for Angled TMA			RP0042

WP Abutments

		Description	GH	Catalog No.
		WP Straight Single TMA Abutment	1.0	WP3731S
			2.0	WP3732S
			3.0	WP3733S
			4.0	WP3729S
			5.0	WP3730S
	WP Retaining Screw for TMA 1mm	WP3406		
	WP Retaining Screw for TMA 2mm	WP3407		
	WP Retaining Screw for TMA 3mm	WP3408		
	WP Retaining Screw for TMA 4mm	WP3409		
	WP Retaining Screw for TMA 5mm	WP3412		

		Description	GH1	GH2	Catalog No.
		WP Angled Single TMA Abutment 17°	1.1	2.5	WP3734S
			2.1	3.5	WP3737S
		WP Angled Single TMA Abutment 30°	1.0	3.5	WP3735S
			1.7	4.1	WP3738S
	WP Retaining Screw for Angled TMA	WP3403			

Transfers Recommended Torque 15(Ncm)

	Description	L	D1	Catalog No.
	Single TMA Impression Coping Open Tray	11.0	4.9	RS5006S
	TMA Screw for Open Tray Transfer			RS5008
	Single TMA Closed Tray Transfer	8.0	4.9	RS5011S
	TMA Screw for Closed Tray Transfer			RS5009

Snap-On Recommended Torque 15(Ncm)

	Description	L	D1	D2	Catalog No.
	Single TMA Snap-On Transfer	8.0	4.9	3.2	RS5012S
	TMA Prosthetic Screw				RS3404
	IRIS™, Impression coping	10	5.5	7.5	IR0179

Analogs Recommended Torque 15(Ncm)

	Description	L	D1	Catalog No.
	Single TMA Abutment Replica	14.6	4.9	RS5004S
	Single TMA Digital Model Replica	9.0	4.9	RS1001S
	Retaining Screw for Digital Model Analog			RS1002

STMA / TMA Scan Bodies Recommended Torque 15(Ncm)

	Description	L	Catalog No.
	Single TMA Scan Body	6.0	TM0006S
		9.0	TM0009S
	TMA Scan Body	6.0	TM0006
		9.0	TM0009
	TMA Scan Body Screw		TM0001

STMA / TMA Healing and Temporary Cylinder Recommended Torque 15(Ncm)



	Description	L	D1	H	Catalog No.
	Single TMA Titanium Temporary	12.0	4.9	10.5	RS4900S
	TMA Titanium Temporary	12.0	4.9	10.5	RS4900
	Single TMA Plastic Casting Sleeve w/Prosthetic Screw	12.0	4.9	10.5	RS5001S
	Single TMA Healing Cap	5.5	4.9	-	RS5005S
	TMA Prosthetic Screw				RS3404

STMA/TMA Cementing Cone Recommended Torque 15(Ncm)




	Description	L	D1	Catalog No.
	Single TMA Cementing Cone, low profile	4.9	5.5	TM0011S
	TMA Cementing Cone, low profile	4.9	5.5	TM0011
	TMA Cementing Cone	4.9	4.9	TM0010
	TMA Scan Body Screw			TM0001

Tool Drivers | For Angled TMA/TMA Prosthetic Screws

For Angled TMA/TMA prosthetic screw

		Description	AH	L	Catalog No.
		Star Hand Driver Short	9.4	23.4	RP0064
		Star Hand Driver Long	14.4	28.4	RP0004
		Star Torque Driver-Short	10.5	20.4	RP0065
		Star Torque Driver-Long	15.5	25.4	RP0005


For straight TMA

		Description	AH	L	Catalog No.
		TMA/FC Hex Torque Driver - Short	10.1	20	RS6197
		TMA/FC Hex Torque Driver - Long	15.1	25	RS6196
		TMA/FC Female Hand Hex Driver-Short	8.4	22.4	RS6193
		TMA/FC Female Hand Hex Driver-Long	14.4	28.4	RS6191
		TMA/FC Hex Handpiece Adapter-Short	8.0	20	RS6194
		TMA/FC Hex Handpiece Adapter-Long	13.0	25	RS6195A



CAD/CAM

Adin's CAD/CAM components and digital libraries enable clinicians and technicians to design and manufacture state-of-the-art restorations supported by Adin's implants, using current digital restoration processes. The variety of components ensures endless restorative options for treatment of single, multiple, partially, or fully edentulous cases.





RP

For Intra-oral and Desktop Scanner	H	Catalog No.	For Medentika PreFace® Holder	H	Ø	Catalog No.
 RP Scan Body	7	RP0088	 RP Ti Blank	15	11.5	RP1211
	10	RP0089		15	15.8	RP1215





WP

For Intra-oral and Desktop Scanner	H	Catalog No.	For Medentika PreFace® Holder	H	Ø	Catalog No.
 WP Scan Body	7	WP0088	 WP Ti Blank	15	11.5	WP1311
	10	WP0089		15	15.8	WP1315

RP Ti Base Abutments



	Description	RH	GH	Catalog No.
	 Ti Base RP Abutment, Engaged	4.0	1.0	TB03411
		4.0	2.0	TB03421
		4.0	3.0	TB03431
		4.0	4.0	TB03441
		6.0	1.0	TB03611
		6.0	2.0	TB03621
		6.0	3.0	TB03631
		6.0	4.0	TB03641
		8.0	1.0	TB03811
		8.0	2.0	TB03821
		8.0	3.0	TB03831
		8.0	4.0	TB03841
	 Ti Base RP Abutment, Non Engaged	4.0	1.0	TB03412
		4.0	2.0	TB03422
		4.0	3.0	TB03432
		4.0	4.0	TB03442
		6.0	1.0	TB03612
		6.0	2.0	TB03622
		6.0	3.0	TB03632
		6.0	4.0	TB03642
		8.0	1.0	TB03812
		8.0	2.0	TB03822
		8.0	3.0	TB03832
		8.0	4.0	TB03842
	 RP Abutment Screw			RP0001

WP Ti Base Abutments

	Description	RH	GH	Catalog No.
	 Ti Base WP Abutment, Engaged	4.0	1.0	TB05411
		4.0	2.0	TB05421
		4.0	3.0	TB05431
		4.0	4.0	TB05441
		6.0	1.0	TB05611
		6.0	2.0	TB05621
		6.0	3.0	TB05631
		6.0	4.0	TB05641
		8.0	1.0	TB05811
		8.0	2.0	TB05821
		8.0	3.0	TB05831
		8.0	4.0	TB05841
	 Ti Base WP Abutment, Engaged Wide	4.0	1.0	TB05411W
		4.0	2.0	TB05421W
		4.0	3.0	TB05431W
		4.0	4.0	TB05441W
		6.0	1.0	TB05611W
		6.0	2.0	TB05621W
		6.0	3.0	TB05631W
		6.0	4.0	TB05641W
		8.0	1.0	TB05811W
		8.0	2.0	TB05821W
		8.0	3.0	TB05831W
		8.0	4.0	TB05841W
	 Ti Base WP Abutment, Non Engaged	4.0	1.0	TB05412
		4.0	2.0	TB05422
		4.0	3.0	TB05432
		4.0	4.0	TB05442
		6.0	1.0	TB05612
		6.0	2.0	TB05622
		6.0	3.0	TB05632
		6.0	4.0	TB05642
		8.0	1.0	TB05812
		8.0	2.0	TB05822
		8.0	3.0	TB05832
		8.0	4.0	TB05842
	 Ti Base WP Abutment, Non Engaged Wide	4.0	1.0	TB05412W
		4.0	2.0	TB05422W
		4.0	3.0	TB05432W
		4.0	4.0	TB05442W
		6.0	1.0	TB05612W
		6.0	2.0	TB05622W
		6.0	3.0	TB05632W
		6.0	4.0	TB05642W
		8.0	1.0	TB05812W
		8.0	2.0	TB05822W
		8.0	3.0	TB05832W
		8.0	4.0	TB05842W
		WP Abutment Screw		WP0001





RP C-Type Ti Base Abutments

For further product information and use please refer to IFU028

	Description	GH	Catalog No.
	 C-Type Ti Base RP Abutment, Engaged	1.0	TB03414
		2.0	TB03424
		3.0	TB03434
		4.0	TB03444
	 C-Type Ti Base RP Abutment, Non Engaged	1.0	TB03415
		2.0	TB03425
		3.0	TB03435
		4.0	TB03445

WP C-Type Ti Base Abutments

For further product information and use please refer to IFU028

	Description	GH	Catalog No.
	 C-Type Ti Base WP Abutment, Engaged	1.0	TB05414
		2.0	TB05424
		3.0	TB05434
		4.0	TB05444
	 C-Type Ti Base WP Abutment, Engaged Wide	1.0	TB05414W
		2.0	TB05424W
		3.0	TB05434W
		4.0	TB05444W
	 C-Type Ti Base WP Abutment, Non Engaged	1.0	TB05415
		2.0	TB05425
		3.0	TB05435
		4.0	TB05445
	 C-Type Ti Base WP Abutment, Non Engaged Wide	1.0	TB05415W
		2.0	TB05425W
		3.0	TB05435W
		4.0	TB05445W




CAD/CAM Angulated Screw Channel

The angulated screw channel (ASC) solution offers a transformative approach to dental restorations. ASC strategically places a screw access hole at an angle within the prosthesis, optimizing both function and esthetics.


By integrating digital planning and specialized drivers, ASC enhances accuracy and predictability in implant-supported restorations. This approach signifies a significant advancement in dental care, prioritizing esthetics and functionality to ensure patient satisfaction.

Adin's angulated screw channel Ti bases support gingival heights of 1-3mm for each implant platform and allow for angulation of up to 20 degrees.


RP Ti Base Abutment Angular

		RH	GH	Catalog No.
	 RP Ti Base Abutment Angular Engaged	4.0	1.5	TB03413
		4.0	2.0	TB03423
		4.0	3.0	TB03433
		6.0	1.5	TB03613
		6.0	2.0	TB03623
		6.0	3.0	TB03633
		8.0	1.5	TB03813
		8.0	2.0	TB03823
		8.0	3.0	TB03833
		RP Abutment Screw		RP0001





WP Ti Base Abutment Angular

		RH	GH	Catalog No.
	 WP Ti Base Abutment Angular Engaged	4.0	1.0	TB05413
		4.0	2.0	TB05423
		4.0	3.0	TB05433
		6.0	1.0	TB05613
		6.0	2.0	TB05623
		6.0	3.0	TB05633
		8.0	1.0	TB05813
		8.0	2.0	TB05823
		8.0	3.0	TB05833
		WP Abutment Screw		WP0001

WP Ti Base Abutment Angular

		RH	GH	Catalog No.
	 WP Ti Base Abutment Angular Engaged	4.0	1.0	TB05413W
		4.0	2.0	TB05423W
		4.0	3.0	TB05433W
		6.0	1.0	TB05613W
		6.0	2.0	TB05623W
		6.0	3.0	TB05633W
		8.0	1.0	TB05813W
		8.0	2.0	TB05823W
		8.0	3.0	TB05833W
		WP Abutment Screw		WP0001

Angulated Tools

		Description	L	AH	Catalog No.
		Angulated Screw Extractor Hand Driver	24.7	10.4	RP6400
		Angulated Star Screw Hand Driver - Short	23.4	9.0	RP6401
		Angulated Star Screw Hand Driver - Long	28.4	14.0	RP6402
		Angulated Star Screw Torque Driver - Short	17.0	11.0	RP6403
		Angulated Star Screw Torque Driver - Long	22.0	16.0	RP6404
		Angulated Star Screw Handpiece Driver - Extra Short	20.0	7.0	RP6405
		Angulated Star Screw Handpiece Driver - Short	25.0	12.0	RP6406
		Angulated Star Screw Handpiece Driver - Long	30.0	17.0	RP6407

TMA™ - Multi Unit Recommended Torque 35(Ncm)

Our TMA™ (Trans Mucosal Abutment) system is indicated for multiple-unit, screw-retained restorations, and may be used in combination with an implant-level framework design.



RP TMA

		Description	GH	Catalog No.
		RP Straight Trans Mucosal Abutment	1.0	RP0032
			2.0	RP0033
			3.0	RP0034
			4.0	RP3729
			5.0	RP3730
		RP Retaining Screw for TMA 1mm		RP0039
		RP Retaining Screw for TMA 2mm		RP0040
		RP Retaining Screw for TMA 3mm		RP0041
		RP Retaining Screw for TMA 4mm		RP3409
		RP Retaining Screw for TMA 5mm		RP3412

		Description	GH1	GH2	Catalog No.
		RP Angled Trans Mucosal Abutment 17°	1.1	2.6	RP3734
		RP Angled Trans Mucosal Abutment 17°	2.3	3.7	RP3735
		RP Angled Trans Mucosal Abutment 30°	1.2	3.7	RP3737
		RP Angled Trans Mucosal Abutment 30°	1.6	4.0	RP3738
		RP Retaining Screw for Angled TMA			RP0042

WP TMA

		Description	GH	Catalog No.
		RP Straight Trans Mucosal Abutment	1.0	WP3731
			2.0	WP3732
			3.0	WP3733
			4.0	WP3729
			5.0	WP3730
		WP Retaining Screw for TMA 1mm		WP3406
		WP Retaining Screw for TMA 2mm		WP3407
		WP Retaining Screw for TMA 3mm		WP3408
		WP Retaining Screw for TMA 4mm		WP3409
		WP Retaining Screw for TMA 5mm		WP3412

		Description	GH1	GH2	Catalog No.
		WP Angled Trans Mucosal Abutment 17° 2.5mm	1.1	2.5	WP3734
		WP Angled Trans Mucosal Abutment 17° 3.5mm	2.1	3.5	WP3737
		WP Angled Trans Mucosal Abutment 30° 3mm	1.0	3.5	WP3735
		WP Angled Trans Mucosal Abutment 30° 4mm	1.7	4.1	WP3738
		WP Angled Trans Mucosal Abutment 52° 5mm	0.9	5.0	WP3739
		WP Angled Trans Mucosal Abutment 60° 5mm	1.4	5.5	WP3740
		WP Retaining Screw for Angled TMA			WP3403

Transfers Recommended Torque 15(Ncm)

	Description	L	D1	Catalog No.
	TMA Impression Coping Open Tray	11.0	4.9	RS5006
	TMA Screw for Open Tray Transfer			RS5008
	TMA Closed Tray Transfer	8.0	4.9	RS5011
	TMA Screw for Closed Tray Transfer			RS5009

Analogs Recommended Torque 15(Ncm)

	Description	L	D1	Catalog No.
	TMA Abutment Replica	16.8	4.9	RS5004
	TMA Digital Model Replica	12.0	4.9	RS1001
	Retaining Screw for Digital Analog			RS1002

TMA Scan Bodies Recommended Torque 15(Ncm)

	Description	L	Catalog No.
	TMA Scan Body	6.0	TM0006
		9.0	TM0009
	TMA Scan Body Screw		TM0001

TMA Healing and Temporary Cylinder Recommended Torque 15(Ncm)

	Description	L	D1	H	Catalog No.
	TMA Titanium Temporary	12.0	4.9	10.5	RS4900
	TMA Plastic Casting Sleeve/Prosthetic Screw	12.0	4.9	10.5	RS5001
	TMA Healing Cap	7.0	4.9	3.5	RS5003
		5.5	4.9	3.5	RS5005
	TMA Prosthetic Screw				RS3404

TMA Prosthetics Recommended Torque 15(Ncm)

	Description	L	D1	D2	Catalog No.
	TMA Cementing Cone, low profile	4.9	5.5	3.2	TM0011
	TMA Cementing Cone	4.9	4.9	3.8	TM0010
	TMA Scan Body Screw				TM0001

Ball Attachment- Overdenture Restoration

Recommended Torque 35(Ncm)



		Description	GH	H	Catalog No.
		RP Ball Attachment	0.5	4.0	RP0051
			1.0	4.5	RP0052
			2.0	5.5	RP0053
			3.0	6.5	RP0054
			4.0	7.5	RP0055
			5.0	8.5	RP0056
		WP Ball Attachment	0.5	4.0	WP0051
			1.0	4.5	WP0052
			2.0	5.5	WP0053
			3.0	6.5	WP0054
			4.0	7.5	WP0055
			5.0	8.5	WP0056

	Description	L	Ø	Catalog No.
	Ball Attachment Replica	16.6	4.0	RS5740

	Description	L	Ø	D2	Catalog No.
	Plastic Ball Cap White (Nylon)	2.6	4.1	3.6	RS2660
	Plastic Ball Cap Pink Soft (Nylon)	2.6	4.1	3.6	RS2662
	Titanium Ball Cap	2.9	4.7	4.7	RS2670
	Stainless Steel Ball Cap	2.90	4.70	4.70	RS2675SS

	Description	Catalog No.
	Torque Driver for Ball Attachment	RS6090


Grip™ - Overdenture Attachment

As an addition to our wide variety of prosthetic solutions.

Implant supported overdentures are among the most practical prosthetic solutions for patients who cannot utilize fixed restorations, or as an interim solution for patients who are waiting for the osseointegration process of their newly placed implants.

The Grip™ attachment system is a versatile solution, designed to ensure maximum denture retention and stability both for maxillary and mandibular dentures. The variety of retentive caps allows dentists to choose the required retention level needed for each individual case and can easily be replaced, when needed.

The Grip™ attachment system is available for Internal Hex in a range of gingival heights, to enable a reliable, long lasting and cost-effective solution for your patients.



Grip™ Abutments

		Description	Catalog No.
		Grip™ RP Abutment 1mm	RP0038
		Grip™ RP Abutment 2mm	RP0057
		Grip™ RP Abutment 3mm	RP0058
		Grip™ RP Abutment 4mm	RP0059
		Grip™ RP Abutment 5mm	RP0060
		Grip™ RP Abutment 6mm	RP0061
		Grip™ WP Abutment 1mm	WP0063
		Grip™ WP Abutment 2mm	WP0057
		Grip™ WP Abutment 3mm	WP0058
		Grip™ WP Abutment 4mm	WP0059
		Grip™ WP Abutment 5mm	WP0060
		Grip™ WP Abutment 6mm	WP0061

	Description	Catalog No.
	Grip™ for TMA	TM0100
	TMA Prosthetic Screw	RS3404
	Grip™ for Bar, M2 thread	RS0101

Included in Sales Pack

Grip™ Processing Set

	Description		Description		Description	Catalog No.
	Housing		Light		Extra Light	RS6426
	Processing		Regular		Block Out Spacer	


Or

Grip™ Extended Range Processing Set

	Description		Description		Description		Description	Catalog No.
	Housing		Light		Zero		Block Out Spacer	RS6427
	Processing		Regular		Extra Light			

Grip™ Retention Caps

Grip Retention 4 Pack			Description	Catalog No.
			Extra Light	RS6417
			Light	RS6418
			Regular	RS6419

Grip™ Extended Range Retention 4 Pack			Description	Catalog No.
			Zero	RS6420
			Extra Light	RS6421
			Light	RS6422
			Regular	RS6423

Grip™ Accessories

	Description	Catalog No.
	Grip™ Impression Copping 4 Pack Includes RS6415 Grip™ Processing 4 Pack	RS6414
	Metal Housing 4 Pack Includes RS6415 Grip™ Processing 4 Pack	RS6413
	Grip™ Block Out Spacer 4 Pack	RS6424
	Grip™ Replica 4 Pack	RS2508
	Grip™ Processing 4 Pack	RS6415
	Grip™ Processing 20 Pack	RS6416

Grip™ Retention Tool

	Description	Catalog No.
	Grip™ core tool	RS6200

Designed for inserting and removing the Grip™ retention males from the metal housing

Surgical Kits

Successful surgical procedures start with having the right instruments and tools to ensure your desired treatment outcomes are met.


That's why our range has everything you need, including drills, keys, ratchets and insertion tools.

Not only do they make procedures safer and simpler, they result in an enhanced patient experience as well.

Our Surgical Kits streamline and simplify your implant procedures.



Premium Kit CF6300




		Drills	Catalog No.			Surgical Tools	Catalog No.
1		Tip Drill	RS8032	18		Surgical Torque Ratchet 35-100 Ncm	RS6111
2		Tri-Step Drill - Long	TD2218	19		Wrench Ratchet Surgical Instruments Handle	RA0050
3		Pilot Drill 2mmD x 8-18mmL	DL1820	20		Depth Probe	RS6180
4		Twist Step Drill 1.8/2.5mmD x 8-18mmL	DL1825				
5		Twist Step Drill 2.3/2.8mmD x 8-18mmL	DL1828			Parallel Pins	Catalog No.
6		Twist Step Drill 2.8/3.2mmD x 8-18mmL	DL1832	21		Parallel Pin Short 10mm (qty. 4)	RS6150
7		Twist Step Drill 3.2/3.6mmD x 8-18mmL	DL1836	22		Parallel Pin Short 16mm (qty. 2)	RS6155
8		Twist Step Drill 3.8/4.2mmD x 8-18mmL	DL1842			Insertion Tools	Catalog No.
9		Twist Step Drill 4.2/4.6mmD x 8-18mmL	DL1846	26		NP 2mm Hex Handpiece Adapter - Short	NP0068
10		Tri-Step Drill - Short	TD1511	27		NP 2mm Hex Handpiece Adapter - Medium	NP0071
11		Pilot Drill 2mmD x 6-13mmL	DL1320	28		NP 2mm Hex Handpiece Adapter - Long	NP0074
12		Twist Step Drill 1.8/2.5mmD x 6-13mmL	DL1325	29		RP 2.25mm Hex Handpiece Adapter - Short	RP0068
13		Twist Step Drill 2.3/2.8mmDx6-13mmL	DL1328	30		RP 2.25mm Hex Handpiece Adapter - Medium	RP0071
14		Twist Step Drill 2.8/3.2mmDx6-13mmL	DL1332	31		RP 2.25mm Hex Handpiece Adapter - Long	RP0074
15		Twist Step Drill 3.2/3.6mmDx6-13mmL	DL1336	32		WP 2.65mm Hex Handpiece Adapter - Short	WP0068
16		Twist Step Drill 3.8/4.2mmDx6-13mmL	DL1342	33		WP 2.65mm Hex Handpiece Adapter - Medium	WP0071
17		Twist Step Drill 4.2/4.6mmDx6-13mmL	DL1346	34		WP 2.65mm Hex Handpiece Adapter - Long	WP0074
						Drill Extension	Catalog No.
				35		Drill Extension	RS9027
						Hex Drivers	Catalog No.
				36		Star Hand Driver - Short	RP0064
				37		Star Hand Driver - Long	RP0004

*Products may be ordered separately as an addition to the kit

*Products may be ordered separately as an addition to the kit

Advanced Kit
CF6302



		Drills	Catalog No.
3		Pilot Drill 2mmD x 8-18mmL	DL1820
4		Twist Step Drill 2.3/2.8mmD x 8-18mmL	DL1828
5		Twist Step Drill 2.8/3.2mmD x 8-18mmL	DL1832
6		Twist Step Drill 3.2/3.6mmD x 8-18mmL	DL1836
7		Twist Step Drill 3.8/4.2mmD x 8-18mmL	DL1842
8		Twist Step Drill 4.2/4.6mmD x 8-18mmL	DL1846

		Surgical Tools	Catalog No.
18		Surgical Torque Ratchet 35-100 Ncm	RS6111
19		Wrench Ratchet Surgical Instruments Handle	RA0050
20		Depth Probe	RS6180

		Parallel Pins	Catalog No.
21		Parallel Pin Short 10mm (qty. 4)	RS6150
22		Parallel Pin Short 16mm (qty. 2)	RS6155


		Insertion Tools	Catalog No.
26		NP 2mm Hex Handpiece Adapter - Short	NP0068
27		NP 2mm Hex Handpiece Adapter - Medium	NP0071
28		NP 2mm Hex Handpiece Adapter - Long	NP0074
29		RP 2.25mm Hex Handpiece Adapter - Short	RP0068
30		RP 2.25mm Hex Handpiece Adapter - Medium	RP0071
31		RP 2.25mm Hex Handpiece Adapter - Long	RP0074
32		WP 2.65mm Hex Handpiece Adapter - Short	WP0068
33		WP 2.65mm Hex Handpiece Adapter - Medium	WP0071
34		WP 2.65mm Hex Handpiece Adapter - Long	WP0074

		Drill Extension	Catalog No.
35		Drill Extension	RS9027

		Hex Drivers	Catalog No.
36		Star Hand Driver - Short	RP0064
37		Star Hand Driver - Long	RP0004

*The advanced kit is also available with short drills | Catalog No. CF6304

Starter Kit CF6303



		Drills	Catalog No.
3		Pilot Drill 2mmD x 8-18mmL	DL1820
4		Twist Step Drill 2.3/2.8mmD x 8-18mmL	DL1828
5		Twist Step Drill 2.8/3.2mmD x 8-18mmL	DL1832
6		Twist Step Drill 3.2/3.6mmD x 8-18mmL	DL1836
7		Twist Step Drill 3.8/4.2mmD x 8-18mmL	DL1842

		Insertion Tools	Catalog No.
26		NP 2mm Hex Handpiece Adapter - Short	NP0068
27		NP 2mm Hex Handpiece Adapter - Medium	NP0071
28		NP 2mm Hex Handpiece Adapter - Long	NP0074
29		RP 2.25mm Hex Handpiece Adapter - Short	RP0068
30		RP 2.25mm Hex Handpiece Adapter - Medium	RP0071
31		RP 2.25mm Hex Handpiece Adapter - Long	RP0074
32		WP 2.65mm Hex Handpiece Adapter - Short	WP0068
33		WP 2.65mm Hex Handpiece Adapter - Medium	WP0071
34		WP 2.65mm Hex Handpiece Adapter - Long	WP0074

		Surgical Tools	Catalog No.
20		Wrench Ratchet	RS6100
21		Wrench Ratchet Surgical Instruments Handle	RA0050

		Hex Drivers	Catalog No.
36		Star Hand Driver - Short	RP0064
37		Star Hand Driver - Long	RP0004

Drills and Tools

Our DLC (Diamond-Like Carbon) Tri-Step/Twist Step Surgical Drills are uniquely designed and coated to improve the surgical drilling process.

Tri-Step and Twist-Step DLC drills simplify the cutting process, save time and minimize risks.

Our improved DLC drill coating provides state-of-the-art hardness and slickness that increases lifespan of the drill, maintains sharpness for extended drilling and generates less heat, thereby reducing chances of heating the bone or complications arising from conventional drilling protocols.

The drills come in the following diameters (mm): 2.0, 2.5, 2.8, 3.2, 3.6, 4.2, 4.6, 5.2, 5.6

Specially designed for Adin's implants and suitable for their simple and convenient drilling protocol.

They are clearly and conveniently marked and can be used with or without drill stops.



Tri-Step drills


Adin's Tri-Step drills condense drilling protocols and avoid the need for using multiple drills throughout complex procedures.

Because friction at the cutting site is reduced, complications common with high speed drilling are decreased and faster and better recovery can occur.


L		Ø	Size	Catalog No.
		1.8-3.2	Short	TD1511
		1.8-3.2	Long	TD2218

Twist Step drills

Short 6-13mm

L		Ø	Color	Catalog No.
		2.0	White	DL1320
		1.8-2.5	Purple	DL1325
		2.3-2.8	Red	DL1328
		2.8-3.2	Blue	DL1332
		3.2-3.6	Green	DL1336
		3.9	Pink	DL1339
		3.8-4.2	Black	DL1342
		4.2-4.6	Yellow	DL1346
		4.8-5.2	Brown	DL1352
		5.2-5.6	Grey	DL1356

Long 8-18mm

L		Ø	Color	Catalog No.
		2.0	White	DL1820
		1.8-2.5	Purple	DL1825
		2.3-2.8	Red	DL1828
		2.8-3.2	Blue	DL1832
		3.2-3.6	Green	DL1836
		3.9	Pink	DL1839
		3.8-4.2	Black	DL1842
		4.2-4.6	Yellow	DL1846
		4.8-5.2	Brown	DL1852
		5.2-5.6	Grey	DL1856

Surgical Tools



	Description	Catalog No.
	Surgical Torque Ratchet 35-100 Ncm	RS6111
	Prosthetic Torque Ratchet 15-35 Ncm	RS6112
	Wrench Ratchet	RS6100
	Surgical Driver	RS7001
	Depth Probe	RS6180

Osteotomes

	Description	Catalog No.
	Osteotome Handle	HA0074
	Osteotome Insert 2.0mmD	HA0061
	Osteotome Insert 2.8mmD	HA0064
	Osteotome Insert 3.25mmD	HA0067
	Osteotome Insert 3.65mmD	HA0070

Tools

	Description	L	AH	Catalog No.
	Wrench Ratchet Surgical Instruments Handle	17	-	RA0050
	Drill Extension	25	11	RS9027
	Parallel Pin	17	10	RS6150
		23	16	RS6155

RP Implant Drivers

	Description	L	AH	Catalog No.
	RP 2.25mm Hex Handpiece Adapter	24.6	5.7	RP0068
		27.6	8.7	RP0071
		36	17.1	RP0074
	RP Handpiece Hex Adapter	20	8	RP0062
		25	13	RP0072
	RP 2.2mm Hex Torque Driver	20	10.1	RP0067
		25	15.1	RP0066

WP Implant Drivers




	Description	L	AH	Catalog No.
	WP 2.65mm Hex Handpiece Adapter	24.8	5.9	WP0068
		27.9	9	WP0071
		36	17.1	WP0074
	WP Handpiece Hex Adapter	20	8	WP0062
		25	13	WP0072
	WP 2.6mm Hex Torque Driver	20	10.1	WP0067
		25	15.1	WP0066

Prosthetic Ratchet




	Description	Catalog No.
	Prosthetic Torque Ratchet 15-35 Ncm	RS6112

Tool Drivers | For Angled TMA/TMA Prosthetic Screws

For Angled TMA/TMA prosthetic screw

		Description	AH	L	Catalog No.
		Star Hand Driver Short	9.4	23.4	RP0064
		Star Hand Driver Long	14.4	28.4	RP0004
		Star Torque Driver-Short	10.5	20.4	RP0065
		Star Torque Driver-Long	15.5	25.4	RP0005
		Star Handpiece Driver - Short	12	30	RP0018
		Star Handpiece Driver - Long	17	35	RP0019

For straight TMA

		Description	AH	L	Catalog No.
		TMA/FC Hex Torque Driver - Short	10.1	20	RS6197
		TMA/FC Hex Torque Driver - Long	15.1	25	RS6196
		TMA/FC Female Hand Hex Driver-Short	8.4	22.4	RS6193
		TMA/FC Female Hand Hex Driver-Long	14.4	28.4	RS6191
		TMA/FC Hex Handpiece Adapter-Short	8.0	20	RS6194
		TMA/FC Hex Handpiece Adapter-Long	13.0	25	RS6195A

	Description	L	AH	Catalog No.
	Torque Driver for Ball Attachment	14.4	4.5	RS6090

Retrieval Screw

	Description	Catalog No.
	RP Retrieval Screw for Abutments	RP3401
	WP Retrieval Screw for Abutments	WP3401

Guided Surgery

The Adinguide™ was created with you in mind, bringing the guided surgery world to your clinic.

The kits were specifically designed to support effortless procedures, accurate predictable implant placement.

The ActiveFlow™ Irrigation Technology is a unique design that forces coolant through the guide to ensure irrigation reaches the bone and prevents overheating.

The Adinguide™ drills are designed with carved channels which allow saline coolant from your handpiece to stream through the sleeve directly into the osteotomy site while drilling.

Our uniquely designed drill kit is intuitive and easy to use by implant diameter and length.

The Tissue Punch and Pilot Drill allow for easy flapless access.

Keyless procedure enables easy handling and full visibility of surgical field.

Self-centering drills controlled by a 5.5mm diameter cylinder at the top of each drill centers the drill in the guide's sleeve.


Self-centering cylinder engages the sleeve before the drill engages the osteotomy in order to ensure correct positioning of the drill.

Built-in stoppers ensure precise and accurate drilling to the desired depth.



Guided Surgery Kit Implants Drills GSK001



	Drills	Catalog No.		Drills	Catalog No.
1		GS5055	5		GS0638
2		GS0625	9		GS0838
3		GS0825	13		GS1038
7		GS1025	17		GS1138
11		GS1125	21		GS1338
15		GS1325	24		GS1638
19		GS1625	6		GS0646
4		GS0632	10		GS0846
8		GS0832	14		GS1046
12		GS1032	18		GS1146
16		GS1132	22		GS1346
20		GS1332	25		GS1646
23		GS1632			



GS Implant Tools
Kit CloseFit™
RP/WP GSK004




Guided surgery template
5.5mmD Sleeve GS0055



CloseFit™ RP/WP
Implant Platform

Guided Surgery Kit CloseFit™ RP/WP Implants Tools GSK004



	Description	Catalog No.		Description	Catalog No.
1	Anchoring Drill S-2.0, 2.0/1.5mmDx19mmL	GS1915	12	Manual Socket L10.5	GS0023
2	Anchoring Drill S-2.0, 2.0/1.7mmDx19mmL	GS1917	13	Manual Socket L14.5	GS0024
3	Handpiece Guide Fixation Screw Driver	GS0028	14	Handpiece Socket L10.5	GS0021
4	Guided RP Handpiece Fixture Mount S-5.5	GS0041	15	Drill Extension	RS9027
5	Guided RP Ratchet Fixture Mount S-5.5	GS0043	16	Star Hand Driver - Long	RP0004
6	Guided RP Ratchet Fixture Mount S-5.5, Screw Retained	GS0037	17	Surgical Torque Ratchet 35-100 Ncm for hex 4mm	GS6111
7	RP Ratchet retrieval tool	GS0046	18	Guide Fixation Screw Driver	GS0027
8	Guided WP Handpiece Fixture Mount S-5.5	GS0042			
9	Guided WP Ratchet Fixture Mount S-5.5	GS0044			
10	Guided WP Ratchet Fixture Mount S-5.5, Screw Retained	GS0039			
11	WP Ratchet retrieval tool	GS0047			



GS Drills Kit
GSK001



Guided surgery template
5.5mmD Sleeve GS0055



CloseFit™ RP/WP
Implant Platform



18mm Drills

	Description	Catalog No.
	Drill S-5.5, 2.5/1.8mmDx18.0mmL	GS1825
	Drill S-5.5, 3.2/2.5mmDx18.0mmL	GS1832
	Drill S-5.5, 3.85/3.2mmDx18.0mmL	GS1838
	Drill S-5.5, 4.6/3.85mmDx18.0mmL	GS1846

Accessories

	Description	Purchase Options	Catalog No.
	Guided Vertical Fixation Screw S-5.5	x3 pack	GS0007
	2.0mmD Anchoring Screw	x3 pack	GS1920
	2.0mmD Anchoring Sleeve	x3 pack	GS0020
	5.5mmD Sleeve	x1 pack	GS0055
		x3 pack	GS0055x3
		x10 pack	GS0055x10

* Illustrations are not to scale

* Availability of the products in some places is limited and depends on local regulations.

* For more information, please contact your local representative.



Industrial Park Alon Tavor, POB 1128, Afula 1811101, Israel
+972-4-6426732 | www.adin-implants.com